

## SafeSheath® Ultra Lite


HEMOSTATIC VALVE PEEL AWAY  
INTRODUCER SYSTEM

# LATEX-FREE

The Adelante® SafeSheath® Ultra Lite Peel Away Introducer System, features a specially engineered valve membrane that allows for reduced insertion forces during procedures.

## Adelante® SafeSheath® Ultra Lite Hemostatic Valve Peel Away Introducer System Features:

- Specially designed membrane allows effortless insertion and removal of ports and central venous catheters
- Consistently clean and smooth peel
- Ergonomically designed, easy splitting breakaway hub
- French size printed on handle and indicated by color-coded introducer cap
- Snap-locking dilator connector
- Decreased risk of blood loss and air embolism
- Standard 13 cm length with a wide range of French sizes available
- Accepts maximum guidewire up to 0.038"/0.96 mm



Model # Kit (Five per box)	Model # Set (Five per box)	French Size	Sheath Length	Dilator Length
ASUP051318	ASUP051318P	5 F	13 cm	18 cm
ASUP061318	ASUP061318P	6 F	13 cm	18 cm
ASUP071318	ASUP071318P	7 F	13 cm	18 cm
ASUP081318	ASUP081318P	8 F	13 cm	18 cm
ASUP851318	ASUP851318P	8.5 F	13 cm	18 cm
ASUP091318	ASUP091318P	9 F	13 cm	18 cm
ASUP101318	ASUP101318P	10 F	13 cm	18 cm
ASUP111318	ASUP111318P	11 F	13 cm	18 cm
N/A	ASUP131318P	13 F	13 cm	18 cm
N/A	ASUP151318P	15 F	13 cm	18 cm
N/A	ASUP161318P	16 F	13 cm	18 cm
N/A	ASUP1651318P	16.5 F	13 cm	18 cm

Introducers supplied sterile. For single use only. Available bulk upon request. Standard lengths are noted above. Other lengths available upon request. For additional information e-mail us at sales@oscor.com

### Adelante® SafeSheath® Ultra Lite Kit Contains:


- One Adelante® SafeSheath® Ultra Lite peel-away sheath
- One dilator
- One guidewire with a flexible straight and J tip with straightener
- One needle with protector
- One syringe

### Adelante® SafeSheath® Ultra Lite Set Contains:

- One Adelante® SafeSheath® Ultra Lite peel-away sheath
- One dilator



\* SafeSheath is a registered trademark of Pressure Products, Inc.



P: 727.937.2511  
F: 727.934.9835  
E: sales@oscor.com  
www.oscor.com

**CAUTION:** Federal Law restricts this device to sale by or on the order of a physician (U.S.A). Please refer to the Instructions for Use for a complete listing of the indications, contraindications, precautions and warnings, where applicable.

