


# Galt|STICK®



Embedded radiopaque  
marker band on all sheaths



Echogenic,  
marked Chiba needle



Echogenic,  
marked Trocar needle

**The GaltStick® Introducer System** is a 6 French Coaxial sheath assembly with stiffening stylet. The GaltStick® allows for the introduction of a .018" guidewire and .038" working guidewire in percutaneous procedures.

- Embedded radiopaque marker band on all sheaths for enhanced visibility
- Smooth dilator/sheath transition
- Variety of kit offerings
- Hydrophilic or silicone coated sheath configurations
- 6Fr outer sheath, 4Fr dilator, 20cm usable length
- Echogenic, marked Trocar needle or Echogenic, marked Chiba needle



## Silicone Coated GaltStick® Introducer Set with Radiopaque Band

| Model Number | Hydrophilic | .018" x 60cm SS/<br>Gold guidewire | .018" x 60cm NT/<br>Gold guidewire | .038" x 145 PTFE 3mm<br>"J" Straight guidewire | Echo<br>Needle | Quantity |
|--------------|-------------|------------------------------------|------------------------------------|------------------------------------------------|----------------|----------|
| KIT-030-04   |             | X                                  |                                    | X                                              | Trocar         | 5        |
| KIT-030-05   |             |                                    | X                                  | X                                              | Trocar         | 5        |
| KIT-030-06   |             | X                                  |                                    |                                                | Trocar         | 10       |
| KIT-030-03   |             |                                    | X                                  |                                                | Trocar         | 10       |

## Hydrophilic Coated GaltStick® Introducer Set with Radiopaque Band

| Model Number | Hydrophilic | .018" x 60cm SS/<br>Gold guidewire | .018" x 60cm NT/<br>Gold guidewire | .038" x 145 PTFE 3mm<br>"J" Straight guidewire | Echo<br>Needle | Quantity |
|--------------|-------------|------------------------------------|------------------------------------|------------------------------------------------|----------------|----------|
| KIT-036-04   | X           | X                                  |                                    | X                                              | Trocar         | 5        |
| KIT-036-05   | X           |                                    | X                                  | X                                              | Trocar         | 5        |
| KIT-036-06   | X           | X                                  |                                    |                                                | Trocar         | 10       |
| KIT-036-03   | X           |                                    | X                                  |                                                | Trocar         | 10       |
| KIT-037-03   | X           |                                    | X                                  |                                                | Chiba          | 10       |
| KIT-037-06   | X           | X                                  |                                    |                                                | Chiba          | 10       |

## Accessories

| Model Number | Description                   | Quantity |
|--------------|-------------------------------|----------|
| NDL-101-00   | 21ga x 15cm Trocar, Echogenic | 5        |
| NDL-109-00   | 22ga x 15cm Chiba, Echogenic  | 5        |
| NDL-109-01   | 22ga x 20cm Chiba, Echogenic  | 5        |

