



Features:

1. *Blue dot leads to precise positioning of TR-closure Band.*
2. *Two different sizes are applicable to almost all patients. Large size is for 20cm to 26cm wrist perimeters while small size is for 16cm to 21cm wrist perimeters.*
3. *Pressure in air balloons is easily adjusted by specially designed syringe tip adaptor.*

Ordering Information:

Ref (Incomplete)	Description
SM-TB-L-TP	Large Size without Syringe, Tyvek-plastic Pouch
SM-TB-S-TP	Small Size without Syringe, Tyvek-plastic Pouch
SM-TB-LS-TP	Large Size with 20ml Syringe, Tyvek-plastic Pouch
SM-TB-SS-TP	Small Size with 20ml Syringe, Tyvek-plastic Pouch

Sunny Medical Device (Shenzhen) Co., Ltd.

4/F and 5/F, Dongpingxing Creative Science Park,
 Longgang District, Shenzhen, Guangdong, P. R. China, 518172
 Tel: +86 755 3383-1052 Fax: +86 755 2263-5989, +852 3171-6516
 Email: sunnymedical_cn@139.com
 Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₀₀₈₆
 ISO13485
 CFDA