



Product Features

Product Benefits

» **Excellent torque control**

Elastic material and integrated design allows 1:1 torque response to deliver the guidewire into the target vessel quickly.

» **Durable and lubricant performance**

Hydrophilic coating offers a durable and smooth approach possible even in tortuous vessels.

» **Preeminent visualization**

Tungsten blended polymer ensures great visualization under X-ray.

Comparison test of lubrication and resistance

Resistance/gf



* Internal bench test data comparison, discrepancies may happens under other standards.

Ordering information

Length/Curve	0.018"/Standard	0.025"/Standard	0.032"/Standard	0.035"/Standard	0.038"/Standard
120cm/A	10A11812	10A12512	10A13212	10A13512	10A13812
120cm/S	10S11812	10S12512	10S13212	10S13512	10S13812
120cm/R1.5mmJ	10J31812	10J32512	10J33212	10J33512	10J33812
120cm/R3mmJ	10J51812	10J52512	10J53212	10J53512	10J53812
150cm/A	10A11815	10A12515	10A13215	10A13515	10A13815
150cm/S	10S11815	10S12515	10S13215	10S13515	10S13815
150cm/R1.5mmJ	10J31815	10J32515	10J33215	10J33515	10J33815
150cm/R3mmJ	10J51815	10J52515	10J53215	10J53515	10J53815
180cm/A	10A11818	10A12518	10A13218	10A13518	10A13818
180cm/S	10S11818	10S12518	10S13218	10S13518	10S13818
180cm/R1.5mmJ	10J31818	10J32518	10J33218	10J33518	10J33818
180cm/R3mmJ	10J51818	10J52518	10J53218	10J53518	10J53818
220cm/A	10A11822	10A12522	10A13222	10A13522	10A13822
220cm/S	10S11822	10S12522	10S13222	10S13522	10S13822
260cm/A	10A11826	10A12526	10A13226	10A13526	10A13826
260cm/S	10S11826	10S12526	10S13226	10S13526	10S13826
260cm/R1.5mmJ	10J31826	10J32526	10J33226	10J33526	10J33826
260cm/R3mmJ	10J51826	10J52526	10J53226	10J53526	10J53826
300cm/A	10A11830	10A12530	10A13230	10A13530	10A13830
300cm/S	10S11830	10S12530	10S13230	10S13530	10S13830
300cm/R1.5mmJ	10J31830	10J32530	10J33230	10J33530	10J33830
300cm/R3mmJ	10J51830	10J52530	10J53230	10J53530	10J53830
	0.025"/Stiff	0.032"/Stiff	0.035"/Stiff	0.038"/Stiff	
150cm/A	10A22515	10A23215	10A23515	10A23815	
150cm/S	—	—	10S23515	10S23815	
180cm/A	10A22518	10A23218	10A23518	10A23818	
180cm/S	—	—	10S23518	10S23818	
220cm/A	—	—	10A23522	10A23822	
260cm/A	—	—	10A23526	10A23826	
260cm/S	—	—	10S23526	10S23826	