

# TORX™

GUIDEWIRE

- 
- Transitionless guidewire
  - Galt Gold® radiopaque tip
  - 1:1 Torque and steerability

## TORX™ Guidewire

Galt Medical Corp. is pleased to introduce the TORX™ Nitinol guidewire. Designed to become your go-to wire, the TORX™ guidewire offers a Nitinol kink-resistant core with a radiopaque Galt Gold® tip for enhanced visibility.



For more information about TORX™ Guidewire or to set up a meeting with a Galt Sales Consultant call (800) 639-2800.

# TORX™

## GUIDEWIRE

### .014" TORX™ Guidewire

| Catalog No. | Diameter | Length | Tip Length | Tip Shape | Tip Angle | Units/Box |
|-------------|----------|--------|------------|-----------|-----------|-----------|
| SGW-038-05  | .014"    | 80     | 5cm        | Angled    | 15°       | 5         |
| SGW-038-06  | .014"    | 180    | 5cm        | Angled    | 15°       | 5         |
| SGW-038-07  | .014"    | 300    | 5cm        | Angled    | 15°       | 5         |

### .018" TORX™ Guidewire

| Catalog No. | Diameter | Length | Tip Length | Tip Shape | Tip Angle | Units/Box |
|-------------|----------|--------|------------|-----------|-----------|-----------|
| SGW-038-00  | .018"    | 80     | 2cm        | Straight  | 0°        | 5         |
| SGW-038-01  | .018"    | 80     | 5cm        | Angled    | 15°       | 5         |
| SGW-038-02  | .018"    | 180    | 5cm        | Angled    | 15°       | 5         |
| SGW-038-08  | .018"    | 60     | 5cm        | Straight  | 0°        | 5         |
| SGW-038-13  | .018"    | 300    | 5cm        | Angled    | 15°       | 5         |

### .025" TORX™ Guidewire

| Catalog No. | Diameter | Length | Tip Length | Tip Shape | Tip Angle | Units/Box |
|-------------|----------|--------|------------|-----------|-----------|-----------|
| SGW-038-03  | .025"    | 180    | 2cm        | Straight  | 0°        | 5         |
| SGW-038-04  | .025"    | 180    | 8cm        | Angled    | 15°       | 5         |
| SGW-038-14  | .025"    | 260    | 8cm        | Angled    | 15°       | 5         |



Transitionless  
guidewire



Galt Gold®  
radiopaque tip



Nitinol core  
1:1 Torque and steerability  
Kink-resistant

