

Arterial Lines

Arterial Catheter Mini-Kits

Arterial Catheter Mini-Kits					
CATALOG NO.	CATHETER	LENGTH	GUIDEWIRE DIAMETER/LENGTH	NEEDLE GAUGE/LENGTH	UNITS
498107	20G	3" (7.6cm)	.021" x 25cm	21G x 1" (2.5cm)	10
498101	20G	6" (15cm)	.021" x 40cm	21G x 2" (5.1cm)	10
498111	18G	3" (7.6cm)	.028" x 25cm	19G x 1.5" (3.8cm)	10
498105	18G	6" (15cm)	.028" x 40cm	19G x 2.25" (5.7cm)	10

Kit Contents

- 1 Radiopaque FEP Catheter With Female Luer Hub
- 1 Guidewire (FC/SS/ST)
- 1 Arterial Needle, Clear Flashback HUB




Arterial Catheter Kits

Arterial Catheter Kits					
CATALOG NO.	CATHETER	LENGTH	GUIDEWIRE DIAMETER/LENGTH	NEEDLE GAUGE/LENGTH	UNITS
498106	20G	3" (7.6cm)	.021" x 25cm	21G x 1" (2.5cm)	10
498100	20G	6" (15cm)	.021" x 40cm	21G x 2" (5.1cm)	10
498110	18G	3" (7.6cm)	.028" x 25cm	19G x 1.5" (3.8cm)	10
498104	18G	6" (15cm)	.028" x 40cm	19G x 2.25" (5.7cm)	10

Kit Contents

- 1 Disposable Towel
- 1 5cc Syringe
- 1 #11 Mini Scalpel
- 3 4" x 4" Gauze Sponges
- 1 9" (23cm) PVC Extension with Stopcock
- 1 Radiopaque FEP Catheter With Female Luer Hub
- 1 000 Silk Suture
- 1 Arterial Needle, Clear Flashback Hub
- 1 Guidewire (FC/SS/ST)



Arterial Catheter Trays

Arterial Catheter Tray					
CATALOG NO.	CATHETER	LENGTH	GUIDEWIRE DIAMETER/LENGTH	NEEDLE GAUGE/LENGTH	UNITS
498411B	20G	3" (7.6cm)	.021" x 25cm	21G x 1" (2.5cm)	10
498410B	20G	6" (15cm)	.021" x 40cm	21G x 2" (5.1cm)	10
498413A	18G	3" (7.6cm)	.028" x 25cm	19G x 1.5" (3.8cm)	10
498412B	18G	6" (15cm)	.028" x 40cm	19G x 2.25" (5.7cm)	10

Kit Contents

- 1 Fenestrated Drape 22" x 22" (56cm x 56cm) With Ripline
- 3 4" x 4" Gauze Sponges
- 4 2" x 2" Gauze Sponges
- 1 Povidone/Iodine Ointment
- 1 000 Silk Suture
- 1 5cc Ampule 1% Lidocaine HCL
- 1 #11 Mini Scalpel
- 1 Radiopaque FEP Catheter With Female Luer Hub
- 1 Arterial Needle, Clear Flashback Hub
- 1 5cc Syringe
- 1 Guidewire (FC/SS/3mm J & ST)
- 1 25G x 5/8" (1.6cm) Needle