

Arterial Lines

- Arterial Catheter Mini-Kits
- Arterial Catheter Kits
- Arterial Catheter Trays
- Arterial Line Kits

Arterial Catheter Mini-Kits

Kit Contents:

- (1) Radiopaque FEP Catheter with Female Luer Hub
- (1) Guidewire (Fixed Core/Stainless Steel/Straight Tip)
- (1) Arterial Needle, Clear Flashback Hub

Box of 10

Catalog #	Catheter	Length	Guidewire		Needle	
			Diameter	Length	Gauge	Length
498107	20G	3" (7.6cm)	.021"	25cm	21G	1" (2.5cm)
498101	20G	6" (15cm)	.021"	40cm	21G	2" (5.1cm)
498111	18G	3" (7.6cm)	.028"	25cm	19G	1.5" (3.8cm)
498105	18G	6" (15cm)	.028"	40cm	19G	2.25" (5.7cm)

Arterial Catheter Kits

Kit Contents:

- (1) Disposable Towel
- (1) 5cc Syringe
- (1) #11 Mini Scalpel
- (3) 4" x 4" Gauze Sponges
- (1) 9" (23cm) PVC Extension with Stopcock
- (1) Radiopaque FEP Catheter with Female Luer Hub
- (1) 000 Silk Suture
- (1) Arterial Needle, Clear Flashback Hub
- (1) Guidewire (Fixed Core/Stainless Steel/Straight Tip)

Box of 10

Catalog #	Catheter	Length	Guidewire		Needle	
			Diameter	Length	Gauge	Length
498106	20G	3" (7.6cm)	.021"	25cm	21G	1" (2.5cm)
498100	20G	6" (15cm)	.021"	40cm	21G	2" (5.1cm)
498110	18G	3" (7.6cm)	.028"	25cm	19G	1.5" (3.8cm)
498104	18G	6" (15cm)	.028"	40cm	19G	2.25" (5.7cm)

Arterial Catheter Trays

Tray Contents:

- (1) Fenestrated Drape 22" x 22" (56cm x 56cm) with Ripline
- (3) 4" x 4" Gauze Sponges
- (4) 2" x 2" Gauze Sponges
- (1) Povidone/Iodine Ointment
- (1) 000 Silk Suture
- (1) 5cc Ampule 1% Lidocaine HCL
- (1) #11 Mini Scalpel
- (1) Radiopaque FEP Catheter with Female Luer Hub
- (1) Arterial Needle, Clear Flashback Hub
- (1) 5cc Syringe
- (1) Guidewire (Fixed Core/Stainless Steel/3mm J and Straight Tip)
- (1) 25G x 5/8" (1.6cm) Needle

Box of 10

Catalog #	Catheter	Length	Guidewire		Needle	
			Diameter	Length	Gauge	Length
498411B	20G	3" (7.6cm)	.021"	25cm	21G	1" (2.5cm)
498410B	20G	6" (15cm)	.021"	40cm	21G	2" (5.1cm)
498412B	18G	6" (15cm)	.028"	40cm	19G	2.25" (5.7cm)

Arterial Line Kits

Kit Contents:

- (1) Disposable Towel
- (1) 5cc Syringe
- (1) #11 Mini Scalpel
- (1) M/F PVC Line LL
- (1) Safety Device
- (3) 4"X4" Gauze Sponges
- (1) 3-way Stopcock
- (1) 000 Silk Suture

Box of 10

Catalog #	Catheter	Length
498443	18G	6" (15cm)
498500	20G	6" (15cm)
498506	20G	3" (7.6cm)