

Drainage

- **SKATER™ Complete Drainage System**
 - SKATER™ All-Purpose and Nephrostomy Drainage Catheters
 - SKATER™ Mini-Loop Drainage Catheters
 - SKATER™ Biliary Drainage Catheters
 - SKATER™ Centesis Catheters
- **SKATER™ Kits**
 - SKATER™ Introducer Sets
- **Accessories**
 - SKATER™ Fix Catheter Fixation Device
 - Drainage Bags
 - Introducer Sheaths/Needles (ISN)
 - Hawkins™ Blunt Needle Access System

SKATER™ Complete Drainage System

The SKATER™ Drainage System provides a complete portfolio of drainage catheters and accessories for all-purpose, nephrostomy, biliary, and centesis procedures. SKATER Drainage Catheters are designed to promote maximum drainage flow and patency, minimize kinking, and are formulated for patient comfort, easy insertion, and radiopacity.

SKATER™ All-Purpose and Nephrostomy Drainage Catheters

SKATER™ All-Purpose and Nephrostomy Drainage Catheters are suitable for percutaneous drainage of abscesses or fluid collections using the direct stick technique or kidneys using the direct stick or Seldinger technique. They can be used safely with ethyl alcohol and doxycycline.


- Each Set Includes:
- (1) SKATER catheter
 - (1) Removable pigtail straightener
 - (1) Metal stiffening cannula
 - (1) Flexible stiffening cannula (not available on 6F and 7F)
 - (1) Locking system with snap-off tab (not available on non-locking pigtail catheters)

SKATER All-Purpose and Nephrostomy Drainage Catheters - Non-Locking Pigtail (Box of 5)

Catalog #	French	Length	Max Guidewire Size	Trocar Stylet
756106015	6F	15cm	.035"	Yes
756106020	6F	20cm	.035"	Yes
756106025	6F	25cm	.035"	Yes
756107020	7F	20cm	.035"	Yes
756107025	7F	25cm	.035"	Yes
756108020	8F	20cm	.038"	Yes
756108025	8F	25cm	.038"	Yes
756108030	8F	30cm	.038"	Yes
756110020	10F	20cm	.038"	Yes
756110025	10F	25cm	.038"	Yes
756110030	10F	30cm	.038"	Yes
756112020	12F	20cm	.038"	Yes
756112025	12F	25cm	.038"	Yes
756114020	14F	20cm	.038"	Yes
756114025	14F	25cm	.038"	Yes
756116025	16F	25cm	.038"	Yes

SKATER All-Purpose and Nephrostomy Drainage Catheters - Locking Pigtail (Box of 5)

Catalog #	French	Length	Max Guidewire Size	Trocar Stylet
756606015	6F	15cm	.035"	Yes
756606020	6F	20cm	.035"	Yes
756606025	6F	25cm	.035"	Yes
756606035	6F	35cm	.035"	No
756607015	7F	15cm	.035"	Yes
756607020	7F	20cm	.035"	Yes
756607025	7F	25cm	.035"	Yes
756608015	8F	15cm	.038"	Yes
756608020	8F	20cm	.038"	Yes
756608025	8F	25cm	.038"	Yes
756608030	8F	30cm	.038"	Yes
756608035	8F	35cm	.038"	No
756608045	8F	45cm	.038"	No
756610015	10F	15cm	.038"	Yes
756610020	10F	20cm	.038"	Yes
756610025	10F	25cm	.038"	Yes
756610030	10F	30cm	.038"	Yes
756610035	10F	35cm	.038"	No
756610045	10F	45cm	.038"	No
756610060	10F	60cm	.038"	No
756612015	12F	15cm	.038"	Yes
756612020	12F	20cm	.038"	Yes
756612025	12F	25cm	.038"	Yes
756612030	12F	30cm	.038"	Yes
756612035	12F	35cm	.038"	No
756612045	12F	45cm	.038"	No
756612060	12F	60cm	.038"	No
756614015	14F	15cm	.038"	Yes
756614020	14F	20cm	.038"	Yes
756614025	14F	25cm	.038"	Yes
756614030	14F	30cm	.038"	Yes
756614035	14F	35cm	.038"	No
756614045	14F	45cm	.038"	No
756614060	14F	60cm	.038"	No
756616025	16F	25cm	.038"	Yes
756616030	16F	30cm	.038"	Yes
756616035	16F	35cm	.038"	No
756616045	16F	45cm	.038"	No
756616060	16F	60cm	.038"	No



SKATER™ Mini-Loop Drainage Catheters


SKATER™ Mini-Loop Drainage Catheters are suitable for percutaneous drainage of abscesses or fluid collections in small cavities using the direct stick or Seldinger technique.

Each Set Includes:

- (1) SKATER catheter
- (1) Removable pigtail straightener
- (1) Metal stylet
- (1) Metal cannula
- (1) Flexible stiffening cannula (not available on 6F and 7F)

SKATER Mini-Loop Drainage Catheters - Locking Pigtail (Box of 5)

Catalog #	French	Length	Max Guidewire Size
758606015	6F	15cm	.035"
758606020	6F	20cm	.035"
758606025	6F	25cm	.035"
758607015	7F	15cm	.035"
758607025	7F	25cm	.035"
758608015	8F	15cm	.038"
758608025	8F	25cm	.038"
758610015	10F	15cm	.038"
758610025	10F	25cm	.038"
758612015	12F	15cm	.038"
758612025	12F	25cm	.038"
758614025	14F	25cm	.038"




SKATER™ Biliary Drainage Catheters

SKATER™ Biliary Drainage Catheters are used to drain bile from the bile duct and are placed using the Seldinger technique.


SKATER Biliary Drainage Catheters - Locking Pigtail (Box of 5)

Catalog #	French	Length	Max Guidewire Size
755308040	8F	40cm	.038"
755310040	10F	40cm	.038"
755312040	12F	40cm	.038"



SKATER™
Centesis Catheters

SKATER™ Centesis Catheters are designed for percutaneous fluid aspirations and small volume drainage procedures.




Box of 5

Catalog #	French	Length
700304007	4F	7cm
700304010	4F	10cm
700304015	4F	15cm
700305007	5F	7cm
700305010	5F	10cm
700305015	5F	15cm
700305020	5F	20cm

SKATER™ Safety Centesis Catheter w/Airlock Valve

SKATER™ Safety Centesis Catheter with Airlock Valve reduces the risk of pneumothorax and fluid leakage.




Box of 10

Catalog #	French	Length
700505010	5F	10cm
700508012	8F	12cm

SKATER™ Safety Centesis Catheter w/Pigtail and Blunt Retracting Stylet

SKATER™ Safety Centesis Catheter with Pigtail and Blunt Retracting Stylet's hemostatic valve minimizes risk of leakage while allowing for reinsertion of a needle or guidewire.




Box of 10

Catalog #	French	Length
700606016	6F	16cm
700608016	8F	16cm

SKATER™ Introducer Sets (ISN)

The SKATER™ Introducer Set provides easy, accurate, and atraumatic placement of a .035"/.038" guidewire in non-vascular procedures.




SKATER Introducer Sets for Introduction of .035" or .038" guidewire (Box of 5)

Catalog #	Description
651506000	4F Dilator/6f X 18cm Coaxial Introducer, .018" x 80cm Stainless Steel Access Guidewire with Radiopaque Tip, 21G x 15cm Trocar Needle
651506300	4FDilator/6f X 18cm Coaxial Introducer, .018" x 80cm Nitinol Access Guidewire with Radiopaque Tip, 21G x 15cm Trocar Needle
651006000	4F Dilator/6f X 18cm Coaxial Introducer, .018" x 80cm Stainless Steel Guidewire with Platinum Tip, 21G x 15cm Trocar Echogenic Needle
651006300	4F Dilator/6f X 18cm Coaxial Introducer, .018" x 80cm Nitinol Guidewire with Platinum Tip, 21G x 15cm Trocar Echogenic Needle

Accessories

SKATER™ Fix Catheter Fixation Device

The SKATER™ Fix Catheter Fixation Device provides suture-free anchoring of percutaneous drainage catheters. It is designed to stop catheter migration and secure the catheter to the patient for up to 7 days.



Box of 25

Catalog #	Description
802500001	SKATER Fix Catheter Fixation Device

Drainage Bags

Designed for efficient fluid collection during drainage procedures.

- Collects up to 600mL of drainage fluid
- Flexible tubing allows patient mobility during drainage
- Twist valve provides quick emptying of drainage bag
- Anti-reflux valve prevents retrograde contamination



Box of 20

Catalog #	Description
DBAG600	600mL Drainage, with Adjustable Drain Tube Length

Box of 5

Catalog #	Description
DBAG600H	600mL Drainage, with Adjustable Drain Tube Length

Introducer Sheaths/Needles (ISN)

The ISN quickly and smoothly provides access for placement of up to a .038" guidewire in non-vascular procedures.

Box of 10

Catalog #	Gauge	Length	Sheath	Max Guidewire Size
ISN1915	19G	15cm	5F	.038"
ISN1920	19G	20cm	5F	.038"

Hawkins™ Blunt Needle Access System

Hawkins™ Blunt Needles offer sharp and blunt stylets for controlled access in percutaneous drainage procedures.

Box of 5

Catalog #	Gauge	Length
216057	16G	5cm
216107	16G	10cm
216157	16G	15cm
218057	18G	5cm
218107	18G	10cm
218157	18G	15cm
218207	18G	20cm

