

Drapes & Covers

- Specialty Drapes
- Angiography Drapes
- Equipment Covers

Specialty Drapes

Our lines of specialty and angiography drapes are made from two materials:

- SMS is a laminated synthetic, nonwoven material that is strong, soft, stable, and repellent; it has the comfort and draping of linen while repelling fluids like a synthetic material
- Polypropylene (poly) is a hydrophobic, strong, semi-soft synthetic material



Box of 20

Catalog #	Description	Material	Size	Absorbent Pad	Fenestration	Transparent Area
397100	EPD Pacing Drape	Poly	75" x 125" (191cm x 318cm)	22" x 125" (56cm x 318cm)	6" x 6.75" (15cm x 17cm)	22" x 125" (56cm x 318cm)
397201	Pacing Drape	SMS	75" x 125" (191cm x 318cm)	32" x 44" (81cm x 112cm)	6" x 6.75" (15cm x 17cm)	
397101	EPD Pacing Drape - Long	Poly	75" x 150" (191cm x 400cm)	22" x 150" (56cm x 400cm)	6" x 6.75" (15cm x 17cm)	22" x 125" (56cm x 318cm)
396260	General Purpose Utility Drape	SMS	40" x 60" (101cm x 152cm)	32" x 34" (81cm x 86cm)	6.5" x 9" (16cm x 23cm)	
396204	Brachial/Axillary/ Cerebral Drape	SMS	38.5" x 60" (98cm x 152cm)	32" x 34" (81cm x 86cm)	3.5" (9cm)	
396302	Minor Procedure Drape with Ripline	SMS	22" x 22" (56cm x 56cm)		3.5" (9cm)	
396303	Heart Biopsy Drape with Facial Port	SMS	77" x 100" (195cm x 254cm)	32" x 34" (81cm x 86cm)	3.5" (9cm)	
396301	Special Procedure	Poly	32" x 34" (81cm x 86cm)	32" x 34" (81cm x 86cm)	3.5" (9cm)	

Angiography Drapes

Box of 20

Catalog #	Description	Material	Size	Absorbent Pad	Fenestration	Transparent Area
396201	Angio Drape PK /1 - Adult	SMS	75" x 125" (191cm x 318cm)	32" x 44" (81cm x 112cm)	5.5" (14cm)	None
396208	Angio Drape PK /8 - Adult	SMS	77" x 110" (195cm x 279cm)	32" x 34" (81cm x 86cm)	5.5" (14cm)	22" x 110" (56cm x 279cm)
396319	Angio Drape PK /19 - Adult	SMS	84" x 125" (213cm x 318cm)	32" x 50" (81cm x 127cm)	5.5" (14cm)	22" x 125" (56cm x 318cm)
396100	EPD Angio Drape - Adult	Poly	75" x 125" (191cm x 318cm)	22" x 125" (56cm x 318cm)	6" (15cm)	22" x 125" (56cm x 318cm)
396110	EPD Basic Pad Drape - Adult	Poly	75" x 125" (191cm x 318cm)	22" x 125" (56cm x 318cm)	6" (15cm)	22" x 125" (56cm x 318cm)
396101	EPD Angio Drape - Neonatal	Poly	75" x 125" (191cm x 318cm)	22" x 125" (56cm x 318cm)	2" (5cm)	22" x 125" (56cm x 318cm)
396102	EPD Angio Drape - Pediatric	Poly	75" x 125" (191cm x 318cm)	22" x 125" (56cm x 318cm)	3" (8cm)	22" x 125" (56cm x 318cm)

Equipment Covers

Box of 20

Catalog #	Description	Width X Height	Features
396223	Bag Type with Elastic Band	20" x 12" (51cm x 30cm)	With Corners
396222	Bag Type with Elastic Band	30" x 17" (76cm x 43cm)	With Corners
396003	Bag Type with Elastic Band	30" x 36" (76cm x 91cm)	With Corners
396221	Bag Type with Elastic Band	36" x 20" (91cm x 51cm)	With Corners
396215	Bag Type with Elastic Band	36" x 28" (91cm x 71cm)	With Corners
396005	Bag Type with Elastic Band	36" x 36" (91cm x 91cm)	With Corners
396007	Bag Type with Elastic Band	42" x 38" (107cm x 96cm)	With Corners
396009	Bag Type with Elastic Band	48" x 36" (122cm x 91cm)	With Corners
396012	Dome Style with Elastic Band	18" Depth (46cm)	
396013	Dome Style with Elastic Band	22" Depth (56cm)	
396397	Dome Style with Elastic Band	27" Depth (69cm)	

Box of 10

Catalog #	Description	Width X Height	Features
396210	Fluoroscope Cover	36" x 44" (91cm x 111cm)	Adhesive on four sides
396348	Back Stand Table Cover	44" x 76" (111cm x 193cm)	