



Endomyocardial Biopsy

- JAWZ™ Endomyocardial Biopsy (EMB) Forceps
- Endomyocardial Biopsy Accessories
- Temporary Pacing Electrode Trays

JAWZ™ Endomyocardial Biopsy (EMB) Forceps

JAWZ™ offers a clean-cut approach to better endomyocardial biopsy samples. It is designed to provide higher quality specimens with fewer passes required. The portfolio consists of disposable devices that maintain sharp dual cutting heads and have great design flexibility.

Dual Cutting Head
Yields Better Samples



- Unique cutting-edge design allows superior surface alignment to provide precise cuts and clean margins
- Precision-machined from a proprietary alloy to remain sharp throughout the procedure

Enhanced Physician Experience

- Designed to deliver higher quality samples with reduced crushed artifact
- Precision-molded to provide a solid "instrument" feel

Greater Flexibility & Ease of Use

- Three-ringed handle allows for full rotations while maintaining your desired thumb position
- Available in Stainless Steel or FEP Coated options and in Straight, Pre-Curved, or Maxi-Curved configurations



Box of 1

Catalog #	Jaw Size	Jaw OD (French)	Jaw Volume	Length	Curve	Coating	Formability	Fenestration
190090	1.5mm	5F	1.80mm ³	50cm	Straight	None	Formable	Fenestrated
190095	1.5mm	5F	1.80mm ³	105cm	Straight	None	Formable	Fenestrated
190051	1.8mm	6F	2.46mm ³	50cm	Maxi-Curved	FEP	Formable	Non-Fenestrated
190060	1.8mm	6F	2.46mm ³	50cm	Straight	None	Formable	Non-Fenestrated
190065	1.8mm	6F	2.46mm ³	50cm	Straight	None	Non- Formable	Non-Fenestrated
190066	1.8mm	6F	2.46mm ³	70cm	Straight	FEP	Formable	Non-Fenestrated
190070	1.8mm	6F	2.46mm ³	105cm	Straight	None	Formable	Non-Fenestrated
190075	1.8mm	6F	2.46mm ³	105cm	Straight	None	Non- Formable	Non-Fenestrated
190080	2.2mm	7F	4.97mm ³	50cm	Pre-Curved	FEP	Formable	Non-Fenestrated
190081	2.2mm	7F	4.97mm ³	50cm	Maxi-Curved	FEP	Formable	Non-Fenestrated
190085	2.2mm	7F	4.97mm ³	105cm	Straight	None	Formable	Non-Fenestrated
190086	2.2mm	7F	4.97mm ³	105cm	Straight	FEP	Formable	Non-Fenestrated
190020	2.4mm	7.5F	5.43mm ³	50cm	Straight	None	Formable	Non-Fenestrated
190030	2.4mm	7.5F	5.43mm ³	50cm	Pre-Curved	FEP	Formable	Non-Fenestrated
190031	2.4mm	7.5F	5.43mm ³	50cm	Maxi-Curved	FEP	Formable	Non-Fenestrated
190040	2.4mm	7.5F	5.43mm ³	105cm	Straight	None	Formable	Non-Fenestrated

Endomyocardial Biopsy Accessories

Heart Biopsy Drape w/Facial Port (Box of 20)

Catalog #	Material	Size	Absorbent Pad	Fenestration
396303	SMS	77" x 100" (195cm x 254cm)	32" x 34" (81cm x 86cm)	3.5" (9cm)

GCI Introducers (Box of 10)

Catalog #	Size	Length	Guidewire Sizes	Hub Color
191545	5F	45cm	.035"-.038"	Gray
191598	5F	98cm	.035"-.038"	Gray
191645	6F	45cm	.035"-.038"	Green
191698	6F	98cm	.035"-.038"	Green
191745	7F	45cm	.035"-.038"	Orange
191798	7F	98cm	.035"-.038"	Orange
191845	8F	45cm	.035"-.038"	Blue
191898	8F	98cm	.035"-.038"	Blue

Temporary Pacing Electrode Trays

Balloon Flow Assisted Bipolar Electrode Tray (Box of 5)

Catalog #	Size	Length
008566A	5F	110cm

Special Care Bipolar Electrode Tray (Box of 5)

Catalog #	Size	Length
008416A	5F	115cm

Notes

[illegible]