


Introducer Sets & Arterial Needles

- **V-Stick™ Vascular Access Sets**
- **Axcess™ Percutaneous Introducer Kits**
- **OneStic™ Micro-Access Introducer Kits**
- **Introducer Kits for Venous Applications**
- **Introducer Kits with Multi-Purpose Curve**
- **Introducer Trays**
- **Specialty Trays & Accessories**
- **Arterial Needles**
- **Guidewire Introducer Needles**
- **Seldinger Needles**
- **Vessel Dilators**

V-Stick™ Vascular Access Sets

V-Stick™ Vascular Access Sets assist in gaining vascular access for the placement of .035" and .038" guidewires into the vascular system using small needle access to reduce puncture site and vessel trauma.

V-Stick Vascular Access Sets include an UltraSharp 21G needle featuring a silicone coated tip and a coaxial introducer, available in standard or stiff, designed to enable smooth transition.



V-Stick Vascular Access Sets w/Standard Coaxial Introducer (Box of 10)

Catalog #	Sheath	Needle			Dilator		Guidewire		
	Size	Gauge	Length	Echogenic	Stiffness	Locking Nut Color	Diameter	Length	Material
391504300E	4F	21G	7cm	Yes	Standard	Red	.018"	40cm	NT
391504400E	4F	21G	7cm	No	Standard	Red	.018"	40cm	NT
391504500E	4F	21G	7cm	No	Standard	Red	.018"	40cm	SS
391505000E	5F	21G	7cm	Yes	Standard	Red	.018"	40cm	SS
391505300E	5F	21G	7cm	Yes	Standard	Red	.018"	40cm	NT
391505400E	5F	21G	7cm	No	Standard	Red	.018"	40cm	NT
391505500E	5F	21G	7cm	No	Standard	Red	.018"	40cm	SS

V-Stick Vascular Access Sets w/Stiff Coaxial Introducer (Box of 10)

Catalog #	Sheath	Needle			Dilator		Guidewire		
	Size	Gauge	Length	Echogenic	Stiffness	Locking Nut Color	Diameter	Length	Material
391604300E	4F	21G	7cm	Yes	Stiff	White	.018"	40cm	NT
391604400E	4F	21G	7cm	No	Stiff	White	.018"	40cm	NT
391604500E	4F	21G	7cm	No	Stiff	White	.018"	40cm	SS
391605000E	5F	21G	7cm	Yes	Stiff	White	.018"	40cm	SS
391605300E	5F	21G	7cm	Yes	Stiff	White	.018"	40cm	NT
391605400E	5F	21G	7cm	No	Stiff	White	.018"	40cm	NT

SS = Stainless Steel
NT = Nitinol

Axcess™ Percutaneous Introducer Kits

Axcess™ Percutaneous Introducer Kits feature a valved introducer to reduce procedural blood loss and air entry. In addition to its ease of use, our Axcess Percutaneous Introducers are available in various lengths and sizes and the kit is available with or without a guidewire.


Axcess Percutaneous Introducer Kit (Standard Length) w/Double-ended Wire (Box of 10)

Catalog #	Sheath		Guidewire		
	Gauge	Length	Hub Color	Diameter	Length
496004	4F	11cm	Red	.035"	45cm
496005	5F	11cm	Gray	.038"	45cm
496006	6F	11cm	Green	.038"	45cm
496007	7F	11cm	Orange	.038"	45cm
496008	8F	11cm	Blue	.038"	45cm
496009	9F	11cm	Black	.038"	45cm



OneStic™ Micro-Access Introducer Kits

OneStic™ Micro-Access Introducer Kits combine Argon’s micro puncture set and valved introducer into one convenient kit. The kit allows for reduced procedure steps with one-step dilation from .018 wire to 4F-8F sheath.



OneStic Micro-Access Introducer Kit w/Stainless Steel Mandrel Wires (Box of 10)

Catalog #	Sheath			Guidewire			Needle		
	Size	Length	Hub Color	Diameter	Length	Core/Tip Material	Gauge	Length	Echogenic
497820	4F	11cm	Red	.018"	45cm	SS/SS	21G	7cm	No
497845	5F	11cm	Gray	.018"	45cm	SS/SS	21G	2.5cm	No
497830	5F	11cm	Gray	.018"	45cm	SS/SS	21G	7cm	No
497840	6F	11cm	Green	.018"	45cm	SS/SS	21G	7cm	No
497861	7F	11cm	Orange	.018"	45cm	SS/SS	21G	7cm	No
497862	8F	11cm	Blue	.018"	45cm	SS/SS	21G	7cm	No

OneStic Micro-Access Introducer Kit w/Nitinol Wires w/Palladium Tip (Box of 10)


Catalog #	Sheath			Guidewire			Needle		
	Size	Length	Hub Color	Diameter	Length	Core/Tip Material	Gauge	Length	Echogenic
497868	4F	5.75cm	Red	.018"	40cm	NT/PD	21G	2.5cm	No
497865	5F	5.75cm	Gray	.018"	40cm	NT/PD	21G	2.5cm	No
497853	6F	5.75cm	Green	.018"	40cm	NT/PD	21G	5.1cm	No
497831	5F	11cm	Gray	.018"	40cm	NT/PD	21G	7cm	No
497841	6F	11cm	Green	.018"	40cm	NT/PD	21G	7cm	No

OneStic Micro-Access Introducer Kit w/Echogenic Needles (Box of 10)

Catalog #	Sheath			Guidewire			Needle		
	Size	Length	Hub Color	Diameter	Length	Core/Tip Material	Gauge	Length	Echogenic
497855	6F	11cm	Green	.018"	40cm	NT/PD	21G	7cm	Yes

SS = Stainless Steel
NT = Nitinol
PD = Palladium

Introducer Kits for Venous Applications



PCI Kit/2 - Introducer Kit w/Detachable Hemostasis Valve (Box of 10)

Catalog #	Sheath		Guidewire	
	Size	Length	Hub Color	Compatibility
496124	8F	11cm	Orange	up to .038"

Kit Contents: (1) Vessel Dilator and (1) Introducer Assembly

PCI Kit/5 - Introducer w/Detachable Hemostasis Valve (Box of 10)

Catalog #	Sheath			Guidewire			
	Size	Length	Hub Color	Diameter	Length	Material	Tip Shade
496183	7F	11cm	Orange	.035"	45cm	SS	Double-ended, 3mm J & Straight

Kit Contents:

- (1) Disposable Towel
- (1) Vessel Dilator
- (3) 4" x 4" Gauze Sponges
- (1) Guidewire (Double-ended, 3mm J, and ST)
- (1) Sheath with Detachable Hemostasis Valve



One-Piece Introducer Kit w/Integral Hemostasis Valve (Box of 10)

Catalog #	Sheath			Guidewire			
	Size	Length	Hub Color	Diameter	Length	Material	Tip Shade
496150	4F	11cm	Red	.025"	45cm	SS	Double-ended, 3mm J & Straight
496151	5F	11cm	Gray	.035"	45cm	SS	Double-ended, 3mm J & Straight
496152	6F	11cm	Green	.035"	45cm	SS	Double-ended, 3mm J & Straight
496153	7F	11cm	Orange	.035"	45cm	SS	Double-ended, 3mm J & Straight
496154	8F	11cm	Blue	.035"	45cm	SS	Double-ended, 3mm J & Straight
496155	8.5F	11cm	White	.035"	45cm	SS	Double-ended, 3mm J & Straight
496156	9F	11cm	Black	.035"	45cm	SS	Double-ended, 3mm J & Straight

SS = Stainless Steel

Introducer Sets
& Arterial Needles

Kit Contents:

- (1) Radiopaque Sheath
- (1) Vessel Dilator
- (1) Guidewire (Double-ended 3mm J and ST)




Tuohy-Borst Introducer Kit w/Integral Hemostasis Valve (Box of 10)

Catalog #	Sheath			Guidewire			
	Size	Length	Hub Color	Diameter	Length	Material	Tip Shade
496132	6F	11cm	Green	.035"	45cm	SS	Double-ended, 3mm J & Straight
496133	7F	11cm	Orange	.035"	45cm	SS	Double-ended, 3mm J & Straight
496134	8F	11cm	Blue	.035"	45cm	SS	Double-ended, 3mm J & Straight
496135	8.5F	11cm	White	.035"	45cm	SS	Double-ended, 3mm J & Straight
496136	9F	11cm	Black	.035"	45cm	SS	Double-ended, 3mm J & Straight

Introducer Kits w/
Multi-Purpose Curve

Box of 10

Catalog #	Sheath			Guidewire
	Size	Length	Hub Color	Compatibility
191545	5F	45cm	Gray	.035" - .038"
191598	5F	98cm	Gray	.035" - .038"
191645	6F	45cm	Green	.035" - .038"
191698	6F	98cm	Green	.035" - .038"
191745	7F	45cm	Orange	.035" - .038"
191798	7F	98cm	Orange	.035" - .038"
191845	8F	45cm	Blue	.035" - .038"
191898	8F	98cm	Blue	.035" - .038"



Introducer Trays

Tray Contents:

- (3) Chloraprep™ Prep Swab
 - (1) Disposable Towel
 - (3) 4" x 4" Gauze Sponges
 - (1) Fenestrated Drape
 - (1) Introducer Sheath
 - (1) 3cc Syringe L/L
 - (1) 22G x 1.5" (3.8cm) Needle
 - (1) #11 Mini Scalpel
 - (1) Septishield Catheter Sleeve
 - (1) 18G x 2.75" (7cm) Arterial Needle, Clear
- (1) Hub
 - (1) 5cc Ampule Lidocaine HCL 1%
 - (1) Guidewire (Double-ended, 3mm J and ST)
 - (2) 5cc Syringe L/S
 - (1) Vessel Dilator
 - (1) 25G x 5/8" (1.6cm) Needle
 - (1) 000 Silk Suture
 - (1) Povidone/Iodine Ointment
 - (1) 18G x 2.5" (6.4cm) Radiopaque O-T-N Catheter



Introducer Tray w/Detachable Hemostasis Valve (Box of 10)


Catalog #	Sheath			Guidewire			
	Size	Length	Hub Color	Diameter	Length	Material	Tip Shade
497422	6F	11cm	Green	.035"	45cm	SS	Double-ended, 3mm J & Straight
497423	7F	11cm	Orange	.035"	45cm	SS	Double-ended, 3mm J & Straight
497424	8F	11cm	Blue	.035"	45cm	SS	Double-ended, 3mm J & Straight
497425	8.5F	11cm	White	.035"	45cm	SS	Double-ended, 3mm J & Straight
497426	9F	11cm	Black	.035"	45cm	SS	Double-ended, 3mm J & Straight

One-Piece Introducer Tray w/Integral Hemostasis Valve (Box of 10)

Catalog #	Sheath			Guidewire			
	Size	Length	Hub Color	Diameter	Length	Material	Tip Shade
497432	6F	11cm	Green	.035"	45cm	SS	Double-ended, 3mm J & Straight
497433	7F	11cm	Orange	.035"	45cm	SS	Double-ended, 3mm J & Straight
497434	8F	11cm	Blue	.035"	45cm	SS	Double-ended, 3mm J & Straight
497436	9F	11cm	Black	.035"	45cm	SS	Double-ended, 3mm J & Straight

SS = Stainless Steel

Specialty Trays & Accessories



Tray Contents:

- (3) Chloraprep™ Prep Swab

(2) Disposable Towel

(6) 4" x 4" Gauze Sponges

(1) Fenestrated Drape w/Ripline

(1) Injection Site

(1) 25G x 5/8" (1.6cm) Needle

(1) 000 Silk Suture

(1) Povidone/Iodine Ointment
- (1) 18G x 2.75" (7cm) Arterial Needle, Clear Hub

(2) 5cc Ampule Lidocaine HCL 1% .035" x 65cm FC, SS, 3mm J & ST

(1) Guidewire

(2) 5cc Syringe L/L

(1) 8F Vessel Dilator

(1) 22G x 1.5" (3.8cm) Needle

(1) #11 Mini Scalpel

(1) 18G x 2.5" (6.4cm) Radiopaque O-T-N Catheter

Universal Percutaneous Introducer Tray (Box of 10)

Catalog #	Description
500200	Universal Percutaneous Introducer Tray

Multi-Purpose Biopsy Kits and Trays (Box of 10)

Catalog #	Description
195100	Multi-Purpose Biopsy Tray

Tray Contents:

- (3) Chloraprep™ Prep Swab

(1) Disposable Towel

(4) 4" x 4" Gauze Sponges

(1) Fenestrated Drape w/Ripline

(1) Plastic Ruler 6" (15cm)

(1) 25G x 5/8" (1.6cm) Needle
- (1) Adhesive Bandage

(6) Microscope Slides

(2) Specimen Tube w/cap

(1) 20G x 1.5" (3.8cm) Needle

(1) 5cc Ampule Lidocaine HCL 1%

(1) 3cc Syringe Luer-Lock
- (1) 5cc Syringe Luer-Slip

(1) 20cc Syringe Luer-Lock

(1) 22G x 1.5" (3.8cm) Needle

(1) #11 Long Handle Scalpel

(1) Medicine Cup 30cc

Accessories (Box of 10)

Catalog #	Description	Size
497296	Detachable Hemostasis Valve with Stopcock	4F - 6F
497297	Detachable Hemostasis Valve with Stopcock	7F - 9F

Accessories (Box of 25)

Catalog #	Description	Size
497302	Septishield Catheter Sleeve (80cm) (Not for use with Axxess™ Introducers)	5F - 8SF



Arterial Needles

AMC/3 - Arterial Needle w/Removable Needle Adapter (Box of 10)

Catalog #	Gauge	Length	Guidewire Compatibility
195312	18G	7cm	up to .038"
195311	21G	2.5cm	up to .021"

AMC/4 - Clear Flashback Needle (Box of 10)


Catalog #	Gauge	Length	Guidewire Compatibility
195317	18G	3.8cm	up to .038"
195318	18G	7cm	up to .038"
195355	19G	7cm	up to .035"
195322	21G	2.5cm	up to .021"
195323	21G	5.1cm	up to .021"

AMC/4 - Winged Arterial Needle (Box of 10)

Catalog #	Gauge	Length	Guidewire Compatibility	Special Features
195352	18G	7cm	up to .038"	
195330	19G	7cm	up to .038"	Extra Thin Wall

Percutaneous Entry Needles (Box of 10)

Catalog #	Description	Gauge	Length	Guidewire Compatibility
195326	Percutaneous Entry Needle (Two-Part Thin Wall Tapered Cannula and Inner Arterial Beveled Stylet)	18G	7.3cm	up to .038"



Guidewire Introducer Needles

Guidewire Introducer Needles are designed for anterior, single wall, arterial percutaneous punctures. The tapered luer hub allows for easy guidewire insertion.

Box of 25

Catalog #	Gauge	Length	Guidewire Compatibility	Special Features
GW11802	18G	7cm	up to .038"	
GW11902	19G	7cm	up to .038"	
GW11802WW	18G	7cm	up to .038"	Clear Window Wall Hub
BGW11802	19G	7cm	up to .038"	Winged Base Plate
BGW11902	19G	7cm	up to .035"	Winged Base Plate

Seldinger Needles

Seldinger Needles are designed for easy, atraumatic, percutaneous arterial punctures.

Box of 25

Catalog #	Gauge	Length	Guidewire Compatibility
SN1802	18G	7cm	up to .038"
SN1902	19G	7cm	up to .038"

Vessel Dilators

Box of 10

Catalog #	Gauge	Length	Guidewire Compatibility	Tip
392305	6F	19.8cm	.035" - .038"	Standard Taper
392306	7F	19.8cm	.035" - .038"	Standard Taper
392307	8F	19.8cm	.035" - .038"	Standard Taper
392309	9F	19.8cm	.035" - .038"	Standard Taper

