

Men's & Women's Health

- HSG Catheters
- Prostate Stabilization Needles
- Prostate Seeding Sets

HSG Catheters

HSG catheters are designed for use in hysterosonography or hysterosalpingography procedures.

Tray Contents:


- (1) HSG Catheter, 5F or 7F
- (1) 18" x 26" Fenestrated Drape
- (1) Welch Allyn® Speculum (compatible with light source #78810)
- (1) Intrauterine Sound
- (2) 3" x 3" Gauze Sponges
- (2) 8" Cotton Swabsticks
- (2) Wide Sponge Swabsticks
- (1) Inner Tray (may be used for Povidone Iodine Solution)
- (1) Povidone Iodine Ointment Packet – 1 gram
- (1) Lubricating Jelly Packet – 5 grams
- (2) Povidone Iodine Solution Packets – 1oz
- (1) Extension Tube Set – 30"
- (1) 20cc Luer-Lock Syringe
- (1) Fixed Volume Syringe
- (1) 18G x 1-1/2" Needle
- (1) Sanitary Napkin
- (1) 30" x 30" CSR Wrap

HSG Catheters (Box of 10)

Catalog #	French	Length	Max Guidewire Size
660005040	5.5F	40cm	.035"
660007040	7F	40cm	.038"

HSG Procedure Trays (Box of 5)

Catalog #	French	Length
660005000	5.5F	40cm
660007000	7F	40cm



Prostate Stabilization Needles

Prostate Stabilization Needles are designed to prevent movement of the prostate gland while either applicator-compatible or preloaded seeding needles are being inserted.



Box of 10 (Sets of 2 needles)

Catalog #	Gauge	Length
500118200	18G	20cm

Prostate Seeding Sets

Prostate Seeding Sets are designed to facilitate transperineal radioactive seed implant procedures.



Box of 25 (4 sets of 5, 5 singles)

Catalog #	Description
500318200SX	For use with Mick® 200 – TP & 200 – TPV applicators 18G x 20cm
500518200X	18G Rigid ETW x 20cm
500618200X	18G Flexible UTW x 20cm