

Introducing the Bengal

Uniquely designed for transradial procedures, the Bengal™ Radial Compression Band offers unprecedented control to achieve radial hemostasis. The simple yet elegant design optimizes targeted compression.

Simple Design

One device in one adjustable size, no syringe or extra gauze required.

Direct Control

Designed to maximize manual control throughout deployment, adjustment and removal

Simple design facilitates easy application and removal, no extra tools required

Transparent polymer provides clear visualization of any bleeding from access site

Fine notches on band allow for micro-adjustment of pressure

Locking mechanism prevents accidental dislodgement during hemostasis period

Thumb-sized space behind the compression pad allows for easy and consistent manual control of the device

Angled structure protects the ulnar nerve and artery from compression

Profile of compression pad delivers a unique pattern: direct pressure on arteriotomy with secondary pressures distal and proximal, designed to prevent hematomas

Adaptor is adjustable for access from the right or left side, and is removable for very small wrists

Proprietary polymer gel ensures patient comfort throughout compression

The Bengal Band includes:

- (1) Compression bracelet
- (1) Adaptor
- (1) Sterile gauze pad

Product Number: RLXX01 Canada and USA


RLXXE01 International (ROW)

Size: Universal

Quantity: 100 per box

Release tab

Ty-rap system





Radial closure device to achieve :

Simple, Quick and effective hemostasis

- No occlusion;
- No discomfort;
- No complications.

