

ACCESSORIES

For vascular and interventional access.

Hemostasis Valve Adapter

Connect with Confidence

- Luer-fit hub, hemostasis valve, & side arm tubing with stopcock
- Can pass up to an 8F device
- Seals on mandrel wires
- Pressure rated up to 30PSI


Gateway)HV®



Needles

Precise Entry and Placement

- Sharp tips for precise access
- X-Sharp 21G echogenic needle features an extra-sharp tip for a smoother, more precise entry
- Bevel-up indicator for visual confirmation of tip orientation



Vessel Dilators

The Pathway for Atraumatic Access


- Atraumatic tips facilitate access and provide a smooth entry with less injury
- Lubricious coating further aids in reducing trauma at the insertion site
- Radiopaque body for visual confirmation of dilator location under fluoroscopy



Guidewire Torque Devices

At Your Command

- Quickly insert, adjust, and remove guidewire with press-and-release finger grips
- Ergonomic, V-shaped design allows for one-handed guidewire manipulation
- Ridges help steer the guidewire tip to the desired location



A Theragenics Company

To find out more, visit
www.theragenics.com/galt



ACCESSORIES

Product Specifications



Percutaneous Entry Thinwall Needles

Catalog No.	Gauge	Length	Needle Tip	Units/Box	
NDL-107-04	21G	21	4cm	Echogenic	25/100
NDL-107-07	21G	21	7cm	Echogenic	25/100
NDL-107-14	21G	21	4cm	Smooth	25/100
NDL-107-17	21G	21	7cm	Smooth	25/100

Percutaneous Entry Thinwall Needles

Catalog No.	Gauge		Length	Needle Tip	Units/Box
NDL-100-00	18G	18	7cm	Smooth	25/100
NDL-100-09	19G	19	7cm	Smooth	25/100
NDL-100-10	19G	19	7cm	Echogenic	25/100
NDL-100-32	21G	21	2.5cm	Smooth	25/100
NDL-100-17	21G	21	4cm	Smooth	25/100
NDL-100-18	21G	21	4cm	Echogenic	25/100
NDL-100-21	21G	21	7cm	Smooth	25/100
NDL-100-22	21G	21	7cm	Echogenic	25/100
NDL-100-28	21G	21	12cm	Echogenic	25/100
NDL-100-29	24G	24	2.5cm	Echogenic	25/100
NDL-100-30	24G	24	4cm	Echogenic	25/100

Vessel Dilators

Catalog No.	Size & Length	Accepts Guidewire Diameter	Units/Box
DIL-600-04	4F 4F x 18cm	.035"	10
DIL-600-05	5F 5F x 18cm	.038"	10
DIL-600-06	6F 6F x 18cm	.038"	10
DIL-600-07	7F 7F x 18cm	.038"	10
DIL-600-08	8F 8F x 18cm	.038"	10
DIL-600-09	9F 9F x 18cm	.038"	10
DIL-600-10	10F 10F x 18cm	.038"	10
DIL-600-11	11F 11F x 18cm	.038"	10
DIL-600-12	12F 12F x 18cm	.038"	10
DIL-600-13	13F 13F x 18cm	.038"	10
DIL-600-14	14F 14F x 18cm	.038"	10
DIL-600-15	15F 15F x 18cm	.038"	10
DIL-600-16	16F 16F x 18cm	.038"	10

QuikTork™

Catalog No.	Accepts Guidewire Diameter	Units/Box
TOR-100-10	.010" - .020" Yellow	10
TOR-100-25	.025" - .040" Orange	10

Gateway)HV*

Catalog No.	Description	Units/Box
GWT-000-30	Hemostasis Valve Adapter	10



A Theragenics Company



galt@theragenics.com



800.639.2800

Theragenics Corporation, 5203 Bristol Industrial Way, Buford, Georgia 30518, USA

Galt Medical Corp. is a subsidiary of Theragenics Corporation

MKT-100-10 Rev. G.