


Designed for Ease of Access

Featuring specially developed angles
and a silicone coating.



Ease of Access



**Up to 57% lower initial force
than the leading competitor***

Fluent transition from cannula, to dilator, to guidewire for
lower penetration forces, improved insertion, and
reduced trauma at the insertion site

Fluent kits available with:


- 4F and 5F sizes with color-coded hubs
- Regular and stiff cannula configurations
- 21G X-Sharp® echogenic needle with bevel indicator
- Variety of guidewire options including high-radiopacity tips
- All-inclusive packaging

4F

5F



Fluent
transition



*1:1 simulated use test data on file at Galt Medical Corp.



A Theragenics Company

To find out more, visit
www.theragenics.com/galt





Fluent Transitionless Micro-Introducer Kits

Catalog No.	Sheath Size		Sheath Length	Guidewire	Guidewire Type	Needle
KIT-038-01	4F	4F	10cm	40cm	SS/Gold	7cm, Echo
KIT-038-02	5F	5F	10cm	40cm	SS/Gold	7cm, Echo
KIT-038-03	4F	4F	10cm	40cm	Nitinol/Gold	7cm, Echo
KIT-038-04	5F	5F	10cm	40cm	Nitinol/Gold	7cm, Echo
KIT-038-05	4F	4F	10cm	60cm	Nitinol/Gold	7cm, Echo
KIT-038-06	5F	5F	10cm	60cm	Nitinol/Gold	7cm, Echo
KIT-038-07	4F	4F	10cm	40cm	SS/SS	7cm, Echo
KIT-038-08	5F	5F	10cm	40cm	SS/SS	7cm, Echo
KIT-038-17	4F	4F	10cm	40cm	Nitinol/Gold	4cm, Echo
KIT-038-18	5F	5F	10cm	40cm	Nitinol/Gold	4cm, Echo
KIT-038-19	4F	4F	10cm	40cm, marked, angled	Nitinol/Gold	7cm, Echo
KIT-038-20	5F	5F	10cm	40cm, marked, angled	Nitinol/Gold	7cm, Echo

Stiffen Fluent Transitionless Micro-Introducer Kits



Catalog No.	Sheath Size		Sheath Length	Guidewire	Guidewire Type	Needle
KIT-039-01	4F	4F	10cm	40cm	SS/Gold	7cm, Echo
KIT-039-03	4F	4F	10cm	40cm	Nitinol/Gold	7cm, Echo
KIT-039-04	5F	5F	10cm	40cm	Nitinol/Gold	7cm, Echo
KIT-039-05	4F	4F	10cm	60cm	SS/SS	7cm, Echo
KIT-039-06	5F	5F	10cm	60cm	SS/SS	7cm, Echo
KIT-039-13	4F	4F	10cm	60cm	Nitinol/Gold	7cm, Echo
KIT-039-14	5F	5F	10cm	60cm	Nitinol/Gold	7cm, Echo
KIT-039-19	4F	4F	10cm	40cm, marked, angled	Nitinol/Gold	7cm, Echo



A Theragenics Company



galt@theragenics.com



800.639.2800

Theragenics Corporation, 5203 Bristol Industrial Way, Buford, Georgia 30518, USA

Galt Medical Corp. is a subsidiary of Theragenics Corporation

MKT-100-07 Rev. I