

Atraumatic Access with Confident Placement

Designed to ease access and confirm placement during drainage, interventional, and other non-vascular percutaneous procedures.


A hydrophilic, coaxial introducer system with dilator and stiffening stylet. Accommodates safety guidewire in parallel with working guidewire.

Atraumatic Access

- Tapered tip and smooth transition minimize injury during insertion
- Lubricious coating further aids in reducing trauma
- Embedded marker band prevents dislodgement for a smooth insertion

Confident Placement

- Stiffening stylet locks in place
- Radiopaque marker band and echogenic needles with depth markers visually confirm tip locations
- Sheath accommodates both .018" safety and .038" working wires
- Gold guidewire tip provides enhanced radiopacity



Stiffening
Stylet

Hydrophilic
coating

Embedded,
radiopaque
Marker band

Tapered tip
with smooth
transition



A Theragenics Company

To find out more, visit
www.theragenics.com/galt



Catalog No.	Sheath/Dilator Size	Guidewire	Echogenic Needle	Box Quantity
KIT-036-02	6F	6F / 4F .018" x 60cm NiTi/Au	N/A	10
KIT-036-04	6F	.018" x 60cm SS/Au .038" x 145cm PTFE 3mm J	Trocar	5
KIT-036-05	6F	.018" x 60cm NiTi/Au .038" x 145cm PTFE 3mm J	Trocar	5
KIT-036-03	6F	.018" x 60cm NiTi/Au	Trocar	10
KIT-036-06	6F	.018" x 60cm SS/Au	Trocar	10
KIT-037-03	6F	.018" x 60cm NiTi/Au	Chiba	10
KIT-037-06	6F	.018" x 60cm SS/Au	Chiba	10

Accessories

Catalog No.	Needle Size	Needle Length	Echogenic Needle	Box Quantity
NDL-101-00	21G	15cm	Trocar	5
NDL-109-00	22G	15cm	Chiba	5
NDL-109-01	22G	20cm	Chiba	5



Our docs love using the GaltStick. We use it for almost every procedure we do in IR!"

Radiologic Technologist



The kit supplies a reliable and consistently accurate wire/needle set. The wire ... is very reliable in following the vessel with ease."

Meghan Darche RT(R)(CT)(VI),
Vascular Institute



The device is nice. Radiopaque band helps with visual. Puncture site is very clean once sheath is removed."

Sarah Postlewaite RT(R),
IR Windy Hill Hospital

Sharp Trocar Needle



Securely locks and aids in penetrating tough anatomy

Sharp Chiba Needle



Allows for aspiration or injection and steers without a stylet



A Theragenics Company



galt@theragenics.com



800.639.2800

Theragenics Corporation, 5203 Bristol Industrial Way, Buford, Georgia 30518, USA

Galt Medical Corp. is a subsidiary of Theragenics Corporation

MKT-100-02 Rev. G