

Micro-Introducers

Micro-access with a Smooth Introduction

Designed for ease of access to facilitate smooth entry into difficult vessels.



A micro-access and Stiffen product line with precise dilator and sheath transitions and lubricious coating.

Ease of Access

- Smooth transition from guidewire to dilator to cannula facilitates vascular access with lower penetration forces and less trauma
- Lubricious coating further aids in reducing trauma at the insertion site
- Entry needle provides tactile feedback and confirms tip location during insertion

Range of Configurations for Procedural Success

- 4F or 5F sheaths in multiple lengths
- 21G smooth or echogenic needles in multiple lengths
- Variety of .018" guidewire options including high-radiopacity tips
- All-inclusive packaging
- Also available in separate sterile components




A Theragenics Company

To find out more, visit
www.theragenics.com/galt



Theragenics

Micro-Introducers



Micro-Introducer Kits

Catalog No.	Sheath Size		Sheath Length	Guidewire	Guidewire Type	Needle	Units/Box
KIT-002-30	4F	4F	10cm	40cm	SS/SS	7cm, Echogenic	10
KIT-002-31	5F	5F	10cm	40cm	SS/SS	7cm, Echogenic	10
KIT-002-36	4F	4F	10cm	60cm	SS/SS	7cm, Echogenic	10
KIT-002-37	5F	5F	10cm	60cm	SS/SS	7cm, Echogenic	10
KIT-002-28	4F	4F	10cm	40cm	SS/Au	7cm, Smooth	10
KIT-002-29	5F	5F	10cm	40cm	SS/Au	7cm, Smooth	10
KIT-002-64	5F	5F	10cm	60cm	NiTi/SS	4cm, Smooth	10
KIT-002-32	4F	4F	10cm	40cm	NiTi/SS	7cm, Echogenic	10
KIT-002-33	5F	5F	10cm	40cm	NiTi/SS	7cm, Echogenic	10
KIT-002-38	4F	4F	10cm	60cm	NiTi/SS	7cm, Echogenic	10
KIT-002-39	5F	5F	10cm	60cm	NiTi/SS	7cm, Echogenic	10
KIT-002-34	4F	4F	10cm	40cm	NiTi/Au	7cm, Echogenic	10
KIT-002-35	5F	5F	10cm	40cm	NiTi/Au	7cm, Echogenic	10
KIT-002-62	4F	4F	10cm	40cm	NiTi/Au	4cm, Smooth	10
KIT-002-63	5F	5F	10cm	40cm	NiTi/Au	4cm, Smooth	10
KIT-002-26	4F	4F	10cm	40cm	SS/Au	7cm, Echogenic	10
KIT-002-27	5F	5F	10cm	40cm	SS/Au	7cm, Echogenic	10

Stiffen Micro-Introducer Kits

Catalog No.	Sheath Size		Sheath Length	Guidewire	Guidewire Type	Needle	Units/Box
KIT-011-02	4F	4F	10cm	40cm	SS/SS	7cm, Echogenic	10
KIT-011-03	5F	5F	10cm	40cm	SS/SS	7cm, Echogenic	10
KIT-011-40	4F	4F	10cm	40cm	NiTi/SS	7cm, Smooth	10
KIT-011-41	5F	5F	10cm	40cm	NiTi/SS	7cm, Smooth	10
KIT-011-42	4F	4F	10cm	40cm	NiTi/SS	7cm, Echogenic	10
KIT-011-43	5F	5F	10cm	40cm	NiTi/SS	7cm, Echogenic	10
KIT-011-62	4F	4F	10cm	40cm	NiTi/Au	7cm, Echogenic	10
KIT-011-63	5F	5F	10cm	40cm	NiTi/Au	7cm, Echogenic	10

Micro-Introducer with Side Holes Kit

Catalog No.	Sheath Size		Sheath Length	Guidewire	Guidewire Type	Needle	Units/Box
KIT-033-04	4F	4F	10cm	40cm	NiTi/SS	7cm, Echogenic	10

Bariatric Micro-Introducer Kit

Catalog No.	Sheath Size		Sheath Length	Guidewire	Guidewire Type	Needle	Units/Box
KIT-011-99	4F	4F	15cm	60cm	NiTi/Au	12cm, Echogenic	10



A Theragenics Company



galt@theragenics.com



800.639.2800

Theragenics Corporation, 5203 Bristol Industrial Way, Buford, Georgia 30518, USA

Galt Medical Corp. is a subsidiary of Theragenics Corporation

MKT-100-17 Rev. F