

Specialty Guidewires

***The precision and design needed
for today's advanced and complex
interventional procedures***

Galt's specialty guidewires are intended for use with diagnostic and therapeutic catheters and interventional devices within the coronary and peripheral vasculature.

Features

- Smooth transition between tip and core
- Stainless, nitinol and PTFE configurations
- Radiopaque gold tip available
- Wide variety of options are available including:
 - Length markings
 - Floppy tips
 - Angled tips



A Theragenics Company

To find out more, visit
www.theragenics.com/galt



Specialty Guidewires

Specialty Guidewires give added performance for today's demanding procedures. Intended for use with diagnostic and therapeutic catheters and interventional devices within the coronary and peripheral vasculature.

Stainless Steel Mandrel – Stainless Steel Coil Guidewire

Catalog No.	Size (Diameter)	Description	Units/Box
SGW-015-02	.018"	40cm	10
SGW-015-08	.018"	45cm	10
SGW-015-03	.018"	60cm	10
SGW-015-14	.018"	80cm	10

Nitinol Mandrel – Stainless Steel Coil Guidewire

Catalog No.	Size (Diameter)	Description	Units/Box
SGW-019-40	.014"	40cm	10
SGW-019-30	.018"	40cm	10
SGW-019-32	.018"	60cm	10
SGW-022-11	.018"	65cm, reference marked	10
SGW-019-34	.018"	80cm	10
SGW-022-13	.018"	130cm, reference marked	10

Nitinol Mandrel – 5.8cm Gold Coil Guidewire

Catalog No.	Size (Diameter)	Description	Units/Box
SGW-056-62	.014"	20cm	5
SGW-056-60	.014"	40cm	5
SGW-056-61	.014"	60cm, angled tip	5
SGW-056-75	.014"	180cm, 4cm angled tip	5
SGW-056-01	.018"	40cm	5
SGW-056-04	.018"	45cm	5
SGW-056-02	.018"	60cm	5
SGW-056-03	.018"	80cm	5

Bentson PTFE Standard Taper Wires

Catalog No.	Size (Diameter)	Description	Units/Box
SGW-035-02	.035"	80cm–20cm taper/6cm distal floppy	10
SGW-035-04	.035"	150cm–20cm taper/6cm distal floppy	10
SGW-035-05	.035"	180cm–20cm taper/6cm distal floppy	10
SGW-035-07	.035"	260cm–20cm taper/6cm distal floppy	10

Stainless Steel PTFE Coated, 3mm “J” Wires

Catalog No.	Size (Diameter)	Description	Units/Box
SGW-011-08	.035"	150cm	10
SGW-011-20	.035"	260cm	10

Note: For .008" – .010" guidewires, see MicroSlide® product line.



A Theragenics Company



galt@theragenics.com



800.639.2800

Theragenics Corporation, 5203 Bristol Industrial Way, Buford, Georgia 30518, USA

Galt Medical Corp. is a subsidiary of Theragenics Corporation

MKT-100-04 Rev. H