

TEARAWAY INTRODUCER SETS & KITS



GALT's Tearaway Introducer Sets & Kits

combine precision ergonomic design with the high quality expected in fine medical devices. Each introducer system provides consistent performance and offers a variety of kit configurations for every need. The smooth transition, easy and secure lock, funnel lead-in and unique patented ergonomic hub design add to a successful access every time.

- 14cm, 30cm, 44cm, and 50cm lengths available
- Sizes available from 4Fr-16Fr, including some half-sizes
- Full kits or individual sets



Tearaway Introducer Sets

Set contents: PTFE Sheath and Dilator only

Catalog Number	Introducer Size	Accepts Wire Diameter	Length	Units/Box
INT-105-02	4.0F	.028"	14cm	10
INT-105-04	5.0F	.038"	14cm	10
INT-105-06	6.0F	.038"	14cm	10
INT-105-07	7.0F	.038"	14cm	10
INT-105-08	8.0F	.038"	14cm	10
INT-105-09	9.0F	.038"	14cm	10
INT-105-11	10.0F	.038"	14cm	10
INT-105-13	11.0F	.038"	14cm	10
INT-105-14	12.0F	.038"	14cm	10
INT-105-16	13.0F	.038"	14cm	10
INT-105-17	14.0F	.038"	14cm	10
INT-105-18	15.0F	.038"	14cm	10
INT-105-19	16.0F	.038"	14cm	10

Long Length Tearaway Introducer Sets

Catalog Number	Introducer Size	Accepts Wire Diameter	Length	Units/Box
INT-106-31	5.0F	.018"	30cm	5
INT-106-09	6.0F	.018"	30cm	5
INT-106-07	7.0F	.018"	30cm	5
INT-106-33	5.0F	.018"	50cm	5
INT-106-37	7.0F	.018"	50cm	5
INT-105-21	5.0F	.038"	30cm	5
INT-105-25	7.0F	.038"	30cm	5
INT-105-29	8.0F	.038"	30cm	5
INT-105-33	9.0F	.038"	30cm	5
INT-105-34	9.0F	.038"	44cm	5
INT-105-37	10.0F	.038"	30cm	5
INT-105-41	11.0F	.038"	30cm	5
INT-105-45	12.0F	.038"	30cm	5
INT-105-53	16.0F	.038"	30cm	5
INT-105-57	6.0F	.038"	30cm	5

5cm, 7cm, and 10cm Micro-Access Introducer Sets also available

Tearaway Introducer Kits

Kit contents: a 14cm Tearaway Introducer Set, 50cm flex tip/3mm "J"-Straight Stainless Guidewire, a 7cm introducer needle, and a 10cc Luer slip syringe

Catalog Number	Introducer Size	Accepts Wire Diameter	Length	Units/Box
KIT-009-03	6.0F	.038"	14cm	5
KIT-009-04	7.0F	.038"	14cm	5
KIT-009-05	8.0F	.038"	14cm	5
KIT-009-07	9.0F	.038"	14cm	5
KIT-009-11	11.0F	.038"	14cm	5

5cm, 7cm, and 10cm Micro-Access Introducer Kits also available



Unique wing design



Easy-to-thread funnel system



Patented locking system



Smooth transition



Galt Medical Corp. • Garland, Texas • Phone: (800) 639-2800 • Fax: (972) 271-4706 • www.galtmedical.com

ISO 13485 Certified. MKT-100-08 Rev. D.