

ISO13485
CE2797
MDR761322 R000
US FDA
NMPA

SUNMED™
昕力®
Since 2003



Sunny Medical

Quality and Affordable
Cardiology and Radiology
Interventional Accessories

- ★ Guide Wire
- ★ Manifold Kit
- ★ Heparin Cap
- ★ Pressure Line
- ★ Control Syringe
- ★ Puncture Needle
- ★ Towel/Sheet Clip
- ★ Introducer Sheath
- ★ Haemostatic Valve
- ★ Angio-Closure Pad
- ★ Manifold and Stopcock
- ★ Pressure Transducer Kit
- ★ Balloon Catheter and SDS
- ★ Balloon Inflation Device (Kit)
- ★ TR-Closure Band/Radiseals

Sunny Medical



公司简介

- 深圳市听力医疗是中国领先的心血管介入手术耗材之生产企业,公司已拥有20多项心血管产品方面的专利。
- 公司拥有三类心脏介入产品生产许可证,并通过了广东省食品药品监督管理局的质量体系考核和获得了欧盟ISO13485质量体系证书。公司13种产品已通过欧盟CE认证—CE2797,并有8种产品通过美国FDA认证—510K。
- 公司的MDD证书已于2023年4月成功转换为MDR证书,听力是中国介入心脏病学界第一个转换MDR成功的公司。截止到2023年4月底,中国仅有不到0.6%的公司成功将MDD证书成功转换为MDR证书。
- 05年公司开始批量出口,目前公司产品已出口到 美国、加拿大、巴西、智利、委内瑞拉、德国、法国、荷兰、西班牙、葡萄牙、瑞士、意大利、波兰、罗马尼亚、保加利亚、塞尔维亚、塞浦路斯、俄罗斯、白俄罗斯、格鲁吉亚、匈牙利、立陶宛、捷克、希腊、土耳其、伊朗、伊拉克、以色列、约旦、黎巴嫩、沙特、也门、阿联酋、南非、埃及、莫桑比克、突尼斯、叙利亚、印度、日本、韩国、巴基斯坦、孟加拉国、泰国、菲律宾、澳大利亚、马来西亚、印尼、新加坡、香港和台湾等50多个国家和地区。
- 公司产品(一次性使用球囊压力泵、一次性使用Y型连接器套件、医用连通器、造影导丝、高压造影注射延长管、桡动脉止血器、股动脉止血器、环柄注射器、穿刺针、肝素帽)已获得中国药监局颁发的证书。

About Us

- Sunny Medical is the leading Contract Manufacturer of medical devices in China, including plastic injection molding, extrusion and assembly, packaging, and ETO sterilization.
- Sunny Medical offers SUNMED™, OEM, ODM and customized medical devices to worldwide customers, specializing in cardiology and radiology intervention products.
- ISO13485 is applied to all the production procedures.
- All mold injection and assembly are in class 100,000 clean room.
- Sterile Test, Positive Reference and Microbial Test are in class 10,000 clean room.
- The MDD certificate is converted to MDR certificate in April 2023, which is the first MDR certificate in China interventional cardiology field. Less than 0.6% of MDD Certificates is successfully converted to MDR based on BSI data till the end of April 2023.
- Sunny has obtained 8 US FDA510K certificates in total. And Sunny is the company with the most FDA510k certificates in the field of cardiovascular interventional in China, including Inflation device, Guide wire, Pressure line, Haemostatic Valves, Control syringes, Disposable Angio-Closure Pads, Disposable needle and Manifold and Stopcocks.
- Two of the top ten medical device companies, M and B, use Sunny's finished products or semi-finished products.

Disposable Balloon Inflation Devices (N-type)



Features:

1. One-Touch release lever for quick deflation with minimal effort.
2. Ergonomic handle designed for easy rotation during inflation and deflation.

Ref (Incomplete)	Description
SM-ID-N3030H-TB	N-type, 30ml, 30atm/bar, HP Stopcock, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +1 949 365 6648
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO13485
FDA510(k) Exempt

Disposable Balloon Inflation Devices (A-type)



Features:

1. Ergonomic design facilitates doctor's operation.
2. Trigger release mechanism provides easy and quick deflation.
3. Inclination of manometer provides easy and precise view of pressure points.


Ref (Incomplete)	Description
SM-ID-A2530H-TB	25ml, 30atm/bar, HP-Stopcock, Tyvek-plastic Blister
SM-ID-A2530HP14-TB	25ml, 30atm/bar, HP Stopcock, with P14 Y connector Kit, Tyvek-plastic Blister
SM-ID-A2530HS14-TB	25ml, 30atm/bar, HP Stopcock, with S14 Y connector Kit, Tyvek-plastic Blister
SM-ID-A2530HC14-TB	25ml, 30atm/bar, HP Stopcock, with C14 Y connector Kit, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE₂₇₉₇
ISO13485
NMPA

Bone Cement Syringe



Features:

1. Spiral design enhancing the product pressure resistance, ensures uniform inflation and smooth injection.
2. Transparent injection cavity with clear scale helps the doctor to observe and give precise injection.
3. Funnel design, being installed in the front end of the cylinder, reduces the loss of bone cement and the risk of contamination.
4. Working pressure: 1500psi.


Ref (Incomplete)	Description
SM-BS-10-TB	10ml, Tyvek-plastic Blister
SM-BS-15-TB	15ml, Tyvek-plastic Blister
SM-BS-20-TB	20ml, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

Disposable Balloon Inflation Devices

Digital display pressure gauge (A-type)



Features:

1. Ergonomic design facilitates doctor's operation.
2. Trigger release mechanism provides easy and quick deflation.
3. Inclination of manometer provides easy and precise view of pressure points.
4. Digital display pressure gauge for precise control of operational pressure.

Ref (Incomplete)	Description
SM-ID-AD2530H-TB	25ml, 30atm/bar(Digital display), HP-Stopcock, Tyvek-plastic Blister
SM-ID-AD2530HP14-TB	25ml, 30atm/bar(Digital display), HP Stopcock, with P14 Y connector Kit, Tyvek-plastic Blister
SM-ID-AD2530HS14-TB	25ml, 30atm/bar(Digital display), HP Stopcock, with S14 Y connector Kit, Tyvek-plastic Blister
SM-ID-AD2530HC14-TB	25ml, 30atm/bar(Digital display), HP Stopcock, with C14 Y connector Kit, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

Disposable Balloon Inflation Devices (F-type)



Features:

1. Accurate inflation and safe deflation.
2. Ergonomic design facilitates doctor's operation.
3. Specially designed automatic cam locking makes it easy-handling.


Ref (Incomplete)	Description
SM-ID-F3030H-TB	F-type, 30ml, 30atm/bar, HP Stopcock, Tyvek-plastic Blister
SM-ID-F3030HP14-TB	F-type, 30ml, 30atm/bar, HP Stopcock, with P14 Y connector Kit, Tyvek-plastic Blister
SM-ID-F3030HS14-TB	F-type, 30ml, 30atm/bar, HP Stopcock, with S14 Y connector Kit, Tyvek-plastic Blister
SM-ID-F3030HC14-TB	F-type, 30ml, 30atm/bar, HP Stopcock, with C14 Y connector Kit, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE₂₇₉₇
ISO13485
FDA510(k) Exempt

Disposable Balloon Inflation Devices (T-type)



Features:

1. Large Fluid Capacity (from 35cc to 60cc).
2. Inflate and deflate as rapidly as 1s.
3. TAVR/TAVI Applicable for valvuloplasty and heart valve replacement.
4. Applicable for esophageal dilation and ERCP procedures.


Ref (Incomplete)	Description
SM-ID-T6020H-TB	T-type, 60ml, 20atm/bar, HP Stopcock, Tyvek-plastic Blister
SM-ID-T5520H-TB	T-type, 55ml, 20atm/bar, HP Stopcock, Tyvek-plastic Blister
SM-ID-T5020H-TB	T-type, 50ml, 20atm/bar, HP Stopcock, Tyvek-plastic Blister
SM-ID-T4520H-TB	T-type, 45ml, 20atm/bar, HP Stopcock, Tyvek-plastic Blister
SM-ID-T4020H-TB	T-type, 40ml, 20atm/bar, HP Stopcock, Tyvek-plastic Blister
SM-ID-T3520H-TB	T-type, 35ml, 20atm/bar, HP Stopcock, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE₂₇₉₇
ISO13485
FDA510(k) Exempt

Haemostatic Valves/Angioplasty Pack (Push-pull Type)



Features:

1. Two different sideons of Y adaptor suit for various customers.
2. Push-pull Y connector makes it easier for doctor to operate precisely.
3. Special valve design to reduce push-pull resistance by 30% plus.
4. Improved valve structure and material reduce moving resistance of wire and catheter by 20% to 50%.

Ref (Incomplete)	Description
SM-HV-7PF-TB	7F, Push-pull, Sideon Female Luer, Tyvek-plastic Blister
SM-HV-8PP-TB	8F, Push-pull, Sideon Tubing and Stopcock, Tyvek-plastic Blister
SM-HV-7PPLCR-TB	7F, Push-pull, Sideon Tubing and Stopcock, Insertion Tool, Red/Copper Torquer, Tyvek-plastic Blister
SM-HV-8PFLCR-TB	8F, Push-pull, Sideon Female Luer, Insertion Tool, Red/Copper Torquer, Tyvek-plastic Blister
SM-HV-7PFLPW-TB	7F, Push-pull, Sideon Female Luer, Insertion Tool, White/Plastic Torquer, Tyvek-plastic Blister
SM-HV-8PPLPW-TB	8F, Push-pull, Sideon Tubing and Stopcock, Insertion Tool, White/Plastic Torquer, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE2797
ISO13485
FDA510(k)
NMPA

Disposable Balloon Inflation Devices (M-type)



Features:

1. Spiral design ensures precise pressure.
2. Improved release handle makes instant and easy inflation and deflation.

Ref (Incomplete)	Description
SM-ID-M2030-TB	M-type, 20ml, 30atm/bar, Tyvek-plastic Blister
SM-ID-M2030H-TB	M-type, 20ml, 30atm/bar, HP Stopcock, Tyvek-plastic Blister
SM-ID-M2030S01-TB	M-type, 20ml, 30atm/bar, with S01 Y connector Kit, Tyvek-plastic Blister
SM-ID-M2030P01-TB	M-type, 20ml, 30atm/bar, with P01 Y connector Kit, Tyvek-plastic Blister
SM-ID-M2030C01-TB	M-type, 20ml, 30atm/bar, with C01 Y connector Kit, Tyvek-plastic Blister
SM-ID-M2030W01-TB	M-type, 20ml, 30atm/bar, with W01 W connector Kit, Tyvek-plastic Blister
SM-ID-M2030PS01-TB	M-type, 20ml, 30atm/bar, with PS01 Y connector Kit, Tyvek-plastic Blister
SM-ID-M2030HS01-TB	M-type, 20ml, 30atm/bar, HP Stopcock, with S01 Y connector Kit, Tyvek-plastic Blister
SM-ID-M2030HP01-TB	M-type, 20ml, 30atm/bar, HP Stopcock, with P01 Y connector Kit, Tyvek-plastic Blister
SM-ID-M2030HC01-TB	M-type, 20ml, 30atm/bar, HP Stopcock, with C01 Y connector Kit, Tyvek-plastic Blister
SM-ID-M2030HW01-TB	M-type, 20ml, 30atm/bar, HP Stopcock, with W01 W connector Kit, Tyvek-plastic Blister
SM-ID-M2030HPS01-TB	M-type, 20ml, 30atm/bar, HP Stopcock, with PS01 Y connector Kit, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE₂₇₉₇

ISO13485

FDA510(k) Exempt
NMPA

Disposable Balloon Inflation Devices (S-type)30ATM/40ATM



Features:

1. Square shape handle applies to different sizes of hands.
2. Chinese knots on handle prevent slippery on holding.
3. Soft cover on deflation release button reduces deflation vibration.


Ref (Incomplete)	Description
SM-ID-S3030H-TB	S-type 30ml, 30atm/bar, HP Stopcock, Tyvek-plastic Blister
SM-ID-S3040H-TB	S-type 30ml, 40atm/bar, HP Stopcock, Tyvek-plastic Blister
SM-ID-S3030HP14-TB	S-type 30ml, 30atm/bar, HP Stopcock, with P14 Y Connector Kit Tyvek-plastic Blister
SM-ID-S3030HS14-TB	S-type 30ml, 30atm/bar, HP Stopcock, with S14 Y Connector Kit, Tyvek-plastic Blister
SM-ID-S3030HC14-TB	S-type 30ml, 30atm/bar, HP Stopcock, with C14 Y Connector Kit, Tyvek-plastic Blister
SM-ID-S3040HP14-TB	S-type 30ml, 40atm/bar, HP Stopcock, with P14 Y Connector Kit Tyvek-plastic Blister
SM-ID-S3040HS14-TB	S-type 30ml, 40atm/bar, HP Stopcock, with S14 Y Connector Kit, Tyvek-plastic Blister
SM-ID-S3040HC14-TB	S-type 30ml, 40atm/bar, HP Stopcock, with C14 Y Connector Kit, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE₂₇₉₇
ISO13485
FDA510(k)

Disposable Balloon Inflation Devices (J-type)



Features:

1. Easy release mechanism ensures one handed operation.
2. Ergonomic handle design provides comfortable grip.
3. Stable pressure gauge with maximum pressure at 40atm.


Ref (Incomplete)	Description
SM-ID-J2830H-TB	28ml, 30atm/bar, HP-Stopcock, Tyvek-plastic Blister
SM-ID-J2840H-TB	28ml, 40atm/bar, HP-Stopcock, Tyvek-plastic Blister
SM-ID-J2830HP14-TB	28ml, 30atm/bar, HP Stopcock, with P14 Y connector Kit, Tyvek-plastic Blister
SM-ID-J2840HP14-TB	28ml, 40atm/bar, HP Stopcock, with P14 Y connector Kit, Tyvek-plastic Blister
SM-ID-J2830HS14-TB	28ml, 30atm/bar, HP Stopcock, with S14 Y connector Kit, Tyvek-plastic Blister
SM-ID-J2840HS14-TB	28ml, 40atm/bar, HP Stopcock, with S14 Y connector Kit, Tyvek-plastic Blister
SM-ID-J2830HC14-TB	28ml, 30atm/bar, HP Stopcock, with C14 Y connector Kit, Tyvek-plastic Blister
SM-ID-J2840HC14-TB	28ml, 40atm/bar, HP Stopcock, with C14 Y connector Kit, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE₂₇₉₇
ISO13485
FDA510(k) Exempt

Haemostatic Valves/Angioplasty Pack (Screw Type)



Features:

1. Two different types of Y adaptor suit for various customers.
2. Large closing knob makes twisting and valve closing very easy.


Ref (Incomplete)	Description
SM-HV-8SF-TB	8F, Screw Type, Sideon Female Luer, Tyvek-plastic Blister
SM-HV-8SP-TB	8F, Screw Type, Sideon Tubing and Stopcock, Tyvek-plastic Blister
SM-HV-8SFLCR-TB	8F, Screw Type, Sideon Female Luer, Insertion Tool, Red/Copper Torquer, Tyvek-plastic Blister
SM-HV-8SPLCR-TB	8F, Screw Type, Sideon Tubing and Stopcock, Insertion Tool, Red/Copper Torquer, Tyvek-plastic Blister
SM-HV-8SFLPW-TB	8F, Screw Type, Sideon Female Luer, Insertion Tool, White/Plastic Torquer, Tyvek-plastic Blister
SM-HV-8SPLPW-TB	8F, Screw Type, Sideon Tubing and Stopcock, Insertion Tool, White/Plastic Torquer, Tyvek-plastic Blister
SM-HV-9SF-TB	9F, Screw Type, Sideon Female Luer, Tyvek-plastic Blister
SM-HV-9SP-TB	9F, Screw Type, Sideon Tubing and Stopcock, Tyvek-plastic Blister
SM-HV-9SFLCR-TB	9F, Screw Type, Sideon Female Luer, Insertion Tool, Red/Copper Torquer, Tyvek-plastic Blister
SM-HV-9SPLCR-TB	9F, Screw Type, Sideon Tubing and Stopcock, Insertion Tool, Red/Copper Torquer, Tyvek-plastic Blister
SM-HV-9SFLPW-TB	9F, Screw Type, Sideon Female Luer, Insertion Tool, White/Plastic Torquer, Tyvek-plastic Blister
SM-HV-9SPLPW-TB	9F, Screw Type, Sideon Tubing and Stopcock, Insertion Tool, White/Plastic Torquer, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE2797
ISO13485
FDA510(k)

Haemostatic Valves/Angioplasty Pack (Y Click)



Features:


1. Three-stage valve creates maneuverability without any bleeding.
2. Auto-closing mechanism ensures the device to move from semi-closed position to closed position automatically when contrast is injected from the side-on lumen.
3. Ergonomic design allows for one hand operation.
4. Improved valve structure and material reduce moving resistance of wire and catheter by 20% to 50%.

Ref (Incomplete)	Description
SM-HV-7CF-TB	7F, Y Click, Sideon Female Luer, Tyvek-plastic Blister
SM-HV-7CP-TB	7F, Y Click, Sideon Tubing and Stopcock, Tyvek-plastic Blister
SM-HV-7CFLCR-TB	7F, Y Click, Sideon Female Luer, Insertion Tool, Red/Copper Torquer, Tyvek-plastic Blister
SM-HV-7CPLCR-TB	7F, Y Click, Sideon Tubing and Stopcock, Insertion Tool, Red/Copper Torquer, Tyvek-plastic Blister
SM-HV-8CFLPW-TB	8F, Y Click, Sideon Female Luer, Insertion Tool, White/Plastic Torquer, Tyvek-plastic Blister
SM-HV-8CPLPW-TB	8F, Y Click, Sideon Tubing and Stopcock, Insertion Tool, White/Plastic Torquer, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

Haemostatic Valves/Angioplasty Pack (Double Screw Type)



Features:

1. Two separate valves in one adaptor fulfill two guide wires operating independently according to actual needs.
2. Different valves make it easier to distinguish and operate the two guide wires.
3. 2-in-1 design saves procedure time and cost.

Ref (Incomplete)	Description
SM-HV-8WF-TB	8F, Double Screw , Sideon Female Luer, Tyvek-plastic Blister
SM-HV-8WFLCR-TB	8F, Double Screw, Sideon Female Luer, Insertion Tool, Red/Copper 2 Torquers, Tyvek-plastic Blister
SM-HV-8WFLPW-TB	8F, Double Screw, Sideon Female Luer, Insertion Tool, White/Plastic 2 Torquers, Tyvek-plastic Blister
SM-HV-9WF-TB	9F, Double Screw, Sideon Female Luer, Tyvek-plastic Blister
SM-HV-9WFLCR-TB	9F, Double Screw, Sideon Female Luer, Insertion Tool, Red/Copper 2 Torquers, Tyvek-plastic Blister
SM-HV-9WFLPW-TB	9F, Double Screw, Sideon Female Luer, Insertion Tool, White/Plastic 2 Torquers, Tyvek-plastic BlisterBlister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE2797
ISO13485
FDA510(k)

Haemostatic Valves/Angioplasty Pack (Short screw Type)



Features:

1. Maximize the availability of catheters.
2. Minimize the interruption by adapters connected.
3. Suitable for different operations.

Ref (Incomplete)	Description
SM-HV-8SSF-TB	8F, Short screw Type, Sideon Female Luer, Tyvek-plastic Blister
SM-HV-8SSL-TB	8F, Short screw Type, Sideon Tubing and Female Luer, Tyvek-plastic Blister
SM-HV-9SSF-TB	9F, Short screw Type, Sideon Female Luer, Tyvek-plastic Blister
SM-HV-9SSL-TB	9F, Short screw Type, Sideon Tubing and Female Luer, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE₂₇₉₇
ISO13485
FDA510(k)

Haemostatic Valves/Angioplasty Pack (DualYst)



Features:

1. Dual-seal hemostasis valve simplifies procedures with an easy one-handed push-and-hold mechanism and allows insertion of interventional devices while controlling bleedback. Valve rated up to 400psi.
2. Dual-seal hemostasis valve simplifies procedures with an easy push-and-hold mechanism, and frees your other hand to deliver interventional devices. Valve rated up to 400psi.


Ref (Incomplete)	Description
SM-HV-7APHF-TB	7F, DualYst, Sideon Female Luer, Sterile, Tyvek-plastic Blister
SM-HV-7APHP-TB	7F, DualYst, Sideon Tubing and Stopcock, Sterile, Tyvek-plastic Blister
SM-HV-7APHFLCR-TB	7F, DualYst, Sideon Female Luer, Insertion Tool with Large Hub, Red/Copper Torquer, Sterile, Tyvek-plastic Blister
SM-HV-7APHPLCR-TB	7F, DualYst, Sideon Tubing and Stopcock, Insertion Tool with Large Hub, Red/Copper Torquer, Sterile, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE₂₇₉₇
ISO13485
FDA510(k)

Haemostatic Valves/Angioplasty Pack (Push-screw Type)



Features:

1. Large 10F inner diameter enables two Balloon Catheters or two *stent* delivery catheters into inner lumen at one time.
2. The use of a hemostatic valve not featured on conventional Y connectors reduces blood loss during procedures.
3. A locking feature of the hemostatic valve has been added, while maintaining the ease-of-use operation.

Ref (Incomplete)	Description
SM-HV-10PSF-TB	10F, Push-screw Type, Sideon Female Luer, Tyvek-plastic Blister
SM-HV-10PSFLCR-TB	10F, Push-screw Type, Sideon Female Luer, Insertion Tool, Red/Copper Torquer, Tyvek-plastic Blister
SM-HV-10PSFLPW-TB	10F, Push-screw Type, Sideon Female Luer, Insertion Tool, White/Plastic Torquer, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE₂₇₉₇
ISO13485
FDA510(k)

Disposable Needles



Features:

1. Sharp tip that has low resistance increases accessibility and reduce tissue trauma.
2. Short tip has better controllability.
3. Tapered hub inner diameter facilitates insertion of the guide wire.
4. Needles of customized sizes are also available.

Ref (Incomplete)	Description
SM-PN-18U07LS-MP	18G, Ultra-Thin Wall, 7cm, Large Hub, Short Tip, Medical Paper-plastic Pouch
SM-PN-18U07SM-MP	18G, Ultra-Thin Wall, 7cm, Small Hub, Medium Tip, Medical Paper-plastic Pouch
SM-PN-18U07BM-MP	18G, Ultra-Thin Wall, 7cm, B Hub, Medium Tip, Medical Paper-plastic Pouch
SM-PN-18T07LM-TB	18G, Thin Wall, 7cm, Large Hub, Medium Tip, Tyvek-plastic Blister
SM-PN-18T07SM-TB	18G, Thin Wall, 7cm, Small Hub, Medium Tip, Tyvek-plastic Blister
SM-PN-18T07BS-TB	18G, Thin Wall, 7cm, B Hub, Short Tip, Tyvek-plastic Blister
SM-PN-21T04LM-MP	21G, Thin Wall, 4cm, Large Hub, Medium Tip, Medical Paper-plastic Pouch
SM-PN-21U04LM-MP	21G, Ultra-Thin Wall, 4cm, Large Hub, Medium Tip, Medical Paper-plastic Pouch
SM-PN-21T04LM-TB	21G, Thin Wall, 4cm, Large Hub, Medium Tip, Tyvek-plastic Blister
SM-PN-21U04LM-TB	21G, Ultra-Thin Wall, 4cm, Large Hub, Medium Tip, Tyvek-plastic Blister
SM-PN-20U05LS-MP	20G, Ultra-Thin Wall, 5cm, Large Hub, Short Tip, Medical Paper-plastic Pouch

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE₂₇₉₇
ISO13485
FDA510(k) Exempt
NMPA

Manifold and Stopcocks



Features:

1. Rotating male luer provides convenient and fast connection.
2. Lubricated rotating handles enhance handle rotatability.


Ref (Incomplete)	Description
SM-MF-25WF2B-TB	250psi, White Handles, Off Handles, 2 Ports, Block Body, Tyvek-plastic Blister
SM-MF-25WO2B-TB	250psi, White Handles, On Handles, 2 Ports, Block Body, Tyvek-plastic Blister
SM-MF-25BO3H-TB	250psi, Blue Handles, On Handles, 3 Ports, Half Body, Tyvek-plastic Blister
SM-MF-25BF3BP-TB	250psi, Blue Handles, Off Handles, 3 Ports, Block Body, Protectors, Tyvek-plastic Blister
SM-MF-50BO3B-TB	500psi, Blue Handles, On Handles, 3 Ports, Block Body, Tyvek-plastic Blister
SM-MF-50BF3B-TB	500psi, Blue Handles, Off Handles, 3 Ports, Block Body, Tyvek-plastic Blister
SM-MF-50WO3B-TB	500psi, White Handles, On Handles, 3 Ports, Block Body, Tyvek-plastic Blister
SM-MF-50WF3B-TB	500psi, White Handles, Off Handles, 3 Ports, Block Body, Tyvek-plastic Blister
SM-MF-50BF3BP-TB	500psi, Blue Handles, Off Handles, 3 Ports, Block Body, Protectors, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE₂₇₉₇
ISO13485
FDA510(k) Exempt
NMPA

Manifold and Stopcocks



Features:

1. Stopcock with rotating male luer can meet high pressure application requirement.


Ref (Incomplete)	Description
SM-SC-04BFSF-MP	400psi, Blue, Off handle, Swivel Nut, Female luer, Medical paper plastic Pouch
SM-SC-05BFRF-TB	500psi, Blue Handle, Off Handle, Rotating Male Luer, Female Luer, Tyvek-Plastic Blister
SM-SC-05WORF-TB	500psi, White Handle, On Handle, Rotating Male Luer, Female Luer, Tyvek-Plastic Blister
SM-SC-12BFRF-TB	1200psi, Blue Handle, Off Handle, Rotating Male Luer, Female Luer, Tyvek-Plastic Blister
SM-SC-12WORF-TB	1200psi, White Handle, On Handle, Rotating Male Luer, Female Luer, Tyvek-Plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE₂₇₉₇
ISO13485
FDA510(k) Exempt
NMPA

Control Syringes



Features:

1. Transparent cylinder enhances visual effect of gauge mark and contrast level.
2. Special designed Silicone piston head and lubrication technology lead plunger easy to move in cylinder and generate high pressure resistance.
3. Thumb and finger rings fullfill precise dose control.

Ref (Incomplete)	Description
SM-CS-10RRR-TB	10ml, Rotating Male Luer, Thumb Ring, Ring Grips, Tyvek-plastic Blister
SM-CS-12RRR-TB	12ml, Rotating Male Luer, Thumb Ring, Ring Grips, Tyvek-plastic Blister
SM-CS-12RRRL-TB	12ml, Rotating Male Luer, Thumb Ring, Ring Grips, Locking Plunger, Tyvek-plastic Blister
SM-CS-10FRR-TB	10ml, Fixed Male Luer, Thumb Ring, Ring Grips, Tyvek-plastic Blister
SM-CS-12FRR-TB	12ml, Fixed Male Luer, Thumb Ring, Ring Grips, Tyvek-plastic Blister
SM-CS-12FRRL-TB	12ml, Fixed Male Luer, Thumb Ring, Ring Grips, Locking Plunger, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE2797
ISO13485
FDA510(k)
NMPA

Control Syringes



Features:

1. Ergonomics plunger design provides comfortable palm pushing.
2. Wing grips design suits for various finger sizes of the doctors and nurses.


Ref (Incomplete)	Description
SM-CS-10RPW-TB	10ml, Rotating Male Luer, Palm Pad, Wing Grips, Tyvek-plastic Blister
SM-CS-12RPW-TB	12ml, Rotating Male Luer, Palm Pad, Wing Grips, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE2797
ISO13485
FDA510(k)
NMPA

Colour Syringes



Features:

1. Multi-color core rod, easy to distinguish the liquid medicine.
2. Precise control of injection dose, simple operation, easy to use.


Ref (Incomplete)	Description
SM-CS-10FR-TB	10ml, Fixed Male Luer, Red, Tyveck-plastic Blister
SM-CS-10FW-TB	10ml, Fixed Male Luer, White, Tyveck-plastic Blister
SM-CS-10FB-TB	10ml, Fixed Male Luer, Blue, Tyveck-plastic Blister
SM-CS-10FY-TB	10ml, Fixed Male Luer, Yellow, Tyveck-plastic Blister
SM-CS-10FG-TB	10ml, Fixed Male Luer, Green, Tyveck-plastic Blister
SM-CS-20FR-TB	20ml, Fixed Male Luer, Red, Tyveck-plastic Blister
SM-CS-20FW-TB	20ml, Fixed Male Luer, White, Tyveck-plastic Blister
SM-CS-20FB-TB	20ml, Fixed Male Luer, Blue, Tyveck-plastic Blister
SM-CS-20FY-TB	20ml, Fixed Male Luer, Yellow, Tyveck-plastic Blister
SM-CS-20FG-TB	20ml, Fixed Male Luer, Green, Tyveck-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE₂₇₉₇
ISO13485
FDA510(k)
NMPA

Disposable Manifold Kits



Features:


1. Saving costs by all-in-one package.
2. Saving time with connected and ready to use kit.

Ref (Incomplete)	Description
SM-MK-FI2C1P1-TP	Manifold with Off and White Handle, Block Body, Control Syringe, 2 Infusion Sets, 1 Pressure Line, Tyvek-plastic Pouch
SM-MK-FI2C1P1-PP	Manifold with Off and White Handle, Block Body, Control Syringe, 2 Infusion Sets, 1 Pressure Line, Plastic-plastic Pouch

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

Disposable TR-closure Band/Radiseals



Features:

1. Blue dot leads to precise positioning of TR-closure Band.
2. Five different sizes are applicable to almost all patients, from XS to XL.
3. One extra extension band can be added to become a universal size.
4. Pressure in air balloons is easily adjusted by specially designed syringe tip adaptor.


Ref (Incomplete)	Description
SM-TB-LS-TP	Large Size with Syringe, Tyvek-plastic Pouch
SM-TB-MS-TP	Medium Size with Syringe, Tyvek-plastic Pouch
SM-TB-SS-TP	Small Size with Syringe, Tyvek-plastic Pouch
SM-TB-XSS-TP	Extra Small Size with Syringe, Tyvek-plastic Pouch
SM-TB-XLS-TP	Extra Large Size with Syringe, Tyvek-plastic Pouch
SM-TB-US-TP	Universal Size with Syringe, Tyvek-plastic Pouch
SM-TB-LBS-TP	Large Size with Syringe, Bonded, Tyvek-plastic Pouch
SM-TB-LLS-TP	Large Size with Syringe, Luer Connector, Tyvek-plastic Pouch
SM-TB-LLB-TP	Large Size without Syringe, Luer Connector, Bonded, Tyvek-plastic Pouch

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE₂₇₉₇
ISO 13485
NMPA

Disposable TR-closure Band/Radiseals



Features:


1. Specially designed 2ml deflation lock mechanism facilitates accurate air discharge.
2. Central balloon design stops bleeding quickly.
3. Bandage on the belt holds the air syringe in place for the doctor to operate.

Ref (Incomplete)	Description
SM-TB-CLLS-TP	C type, Large Size with Syringe, Luer Connector, Tyvek-plastic Pouch
SM-TB-CMLS-TP	C type, Medium Size with Syringe, Luer Connector , Tyvek-plastic Pouch
SM-TB-CSLS-TP	C type, Small Size with Syringe, Luer Connector , Tyvek-plastic Pouch

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

Guide Wires



Features:

1. Flexible tip makes tip easy to deflect.
2. Elastic and shape-resistance distal core and lubricious coating provide great accessibility and excellent steerability.
3. Supportive proximal shaft increases pushability.


Ref (Incomplete)	Description
SM-GW-3S32260P-TP	S-tip, 0.032"x260cm, PTFE Coated, with M-type dispenser, Tyvek-plastic Pouch
SM-GW-3J32260P-TP	3mm Radius, J-tip, 0.032"x260cm, PTFE Coated, with M-type dispenser, Tyvek-plastic Pouch
SM-GW-3S35150P-TP	S-tip, 0.035"x150cm, PTFE Coated, with M-type dispenser, Tyvek-plastic Pouch
SM-GW-3J35150P-TP	3mm Radius, J-tip, 0.035"x150cm, PTFE Coated, with M-type dispenser, Tyvek-plastic Pouch
SM-GW-3S35180P-TP	S-tip, 0.035"x180cm, PTFE Coated, with M-type dispenser, Tyvek-plastic Pouch
SM-GW-3J35180P-TP	3mm Radius J-tip, 0.035"x180cm, PTFE Coated, with M-type dispenser, Tyvek-plastic Pouch
SM-GW-3D35180P-TP	3mm Radius D-tip, 0.035"x180cm, PTFE Coated, with M-type dispenser, Tyvek-plastic Pouch
SM-GW-3S35260PN-TP	S-tip, 0.035"x260cm, PTFE Coated, with N-type dispenser, Tyvek-plastic Pouch
SM-GW-3J35260PN-TP	3mm Radius, J-tip, 0.035"x260cm, PTFE Coated, with N-type dispenser, Tyvek-plastic Pouch
SM-GW-3J38260PN-TP	3mm Radius, J-tip, 0.038"x260cm, PTFE Coated, with N-type dispenser, Tyvek-plastic Pouch
SM-GW-3A35150H-TP	3cm flexible length A-tip, 0.035"x150cm, Hydrophilic Coated, with N-type dispenser, Tyvek-plastic Pouch
SM-GW-3A35180H-TP	3cm flexible length A-tip, 0.035"x180cm, Hydrophilic Coated, with N-type dispenser, Tyvek-plastic Pouch
SM-GW-3A35260H-TP	3cm flexible length A-tip, 0.035"x260cm, Hydrophilic Coated, with N-type dispenser, Tyvek-plastic Pouch
SM-GW-3J14190HF-TP	3cm Radiopacity J-tip, 0.014"x190cm, Hydrophilic Coated, Floppy Tip and Standard Shaft, with N-type dispenser, Tyvek-plastic Pouch
SM-GW-3J14300HF-TP	3cm Radiopacity J-tip, 0.014"x300cm, Hydrophilic Coated, Floppy Tip and Standard Shaft, with N-type dispenser, Tyvek-plastic Pouch

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE2797
ISO13485
FDA510(k)
NMPA

Guide Wires Accessories



Guide wire dispenser straightener

Guide wire dispenser adaptor

Guide wire without coating




Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE₂₇₉₇
ISO13485
FDA510(k)

Disposable Angio-closure Pads



Features:

1. Luer lock fits both ordinary syringe and control syringe.
2. Tape with super adhesive guarantees air balloon having sufficient stanching pressure.


Ref (Incomplete)	Description
SM-AP-L-MP	Large size, 30cmx27.5cm, Medical paper-plastic pouch
SM-AP-L150-MP	Large size, 150cc, Medical paper-plastic pouch
SM-AP-L150S-MP	Large size with Syringe, 150cc, Medical paper-plastic pouch

Sunny Medical Device (Shenzhen) Co., Ltd.


1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE₂₇₉₇
ISO13485
FDA510(k)
NMPA

Banded Bag with Rubber Band and Tape



Banded Bag with Elastic Opening, Circular




Ref (Incomplete)	Description
SM-EB-R91102	Banded bag with rubber band and Tape, 91cm (36in)*102cm(40in)
SM-EB-C50	Banded bag with elastic opening, circular, 50cm (20in)*50cm(20in)
SM-EB-C102	Banded bag with elastic opening, circular, 102cm (40in)*102cm(40in)
SM-EB-C127	Banded bag with elastic opening, circular, 127cm (50in)*127cm(50in)

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

Pressure Lines



Features:

1. 400psi/500psi/900psi/1200psi pressure line by PVC Tubing.
2. Low cost, high performance, without DEHP.
3. 1200psi Pressure Line by Braided PU Tubing.
4. Flexible tubing and rotating male luer provide easy and fast connection.

Ref (Incomplete)	Description
SM-PL-04PC120MFP-MP	400psi, PVC, Coiled,120cm, Male to Female, Protectors,Medical paper-plastic Pouch
SM-PL-04PC120MF-MP	400psi, PVC, Coiled,120cm, Male to Female, Medical paper-plastic Pouch
SM-PL-05P120MF-MP	500psi, PVC, 120cm, Male to Female, Medical paper-plastic PouchBlister
SM-PL-12P120MF-MP	1200psi, PVC, 120cm, Male to Female, Medical paper-plastic Pouch
SM-PL-12P120RF-MP	1200psi, PVC, 120cm, Rotating Male to Female, Medical paper-plastic Pouch
SM-PL-12B120RF-TP	1200psi, Braided PU, 120cm, Rotating Male to Female, Tyvek-plastic Pouch

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE2797
ISO13485
FDA510(k)
NMPA

Accessories



Rotating male luer



Fixed male luer



Fixed female luer



Small hub



Large hub



B hub



Large hub with wing

Clipless Dispenser



Stopper



Protector



Sunny Medical Device (Shenzhen) Co., Ltd.


1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172

Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516

Email: sunnymedical_cn@139.com jamesqizhang@gmail.com

Website: www.sunnymedical.com.cn www.sunnymedical.com

Disposable Invasive Blood Pressure Transducer



Features:


1. Superior visualization of fluid.
2. Consistent and accurate monitoring.
3. Easy assembly, reduce setup time.
4. High quality chip from global leading sensor supplier, Measurement Specialities.

Ref (Incomplete)	Description
SM-DT-UT1-TB	1 Transducer with Utah Type of Cable Connector (TCC), Tyvek-plastic Blister
SM-DT-BB1-TB	1 Transducer with B Braun TCC, Tyvek-plastic Blister
SM-DT-BD1-TB	1 Transducer with BD TCC, Tyvek-plastic Blister
SM-DT-ED1-TB	1 Transducer with Edward TCC, Tyvek-plastic Blister
SM-DT-AB1-TB	1 Transducer with Abbott TCC, Tyvek-plastic Blister
SM-DT-UT1AP2S1-TB	1 Transducer with Utah TCC, 1 Administration Set, 2 Pressure Lines, One Stopcock, Tyvek-plastic Blister
SM-DT-BB1AP2S1-TB	1 Transducer with B Braun TCC, 1 Administration Set, 2 Pressure Lines, One Stopcock, Tyvek-plastic Blister
SM-DT-BD1AP2S1-TB	1 Transducer with BD TCC, 1 Administration Set, 2 Pressure Lines, One Stopcock, Tyvek-plastic Blister
SM-DT-ED1AP2S1-TB	1 Transducer with Edward TCC, 1 Administration Set, 2 Pressure Lines, One Stopcock, Tyvek-plastic Blister
SM-DT-AB1AP2S1-TB	1 Transducer with Abbott TCC, 1 Administration Set, 2 Pressure Lines, One Stopcock, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

Disposable Invasive Blood Pressure Transducer




Ref (Incomplete)	Description
OEM-DT-P-NBP	Plate for Transducer
OEM-DT-C-NBP	Clamp for Transducer
OEM-DT-AC-UT2MR	IBP Adapter Cable Utah Pressure Transducer to Mindray monitor.
OEM-DT-AC-UT2HP	IBP Adapter Cable Utah Pressure Transducer to Philips/HP monitor.
OEM-DT-AC-AB2MR	IBP Adapter Cable Abbott Pressure Transducer to Mindray monitor.
OEM-DT-AC-AB2HP	IBP Adapter Cable Abbott Pressure Transducer to Philips/HP monitor.

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

Balloon Dilatation Catheter



- *Ideal for CTO*
- *Semi-Compliance*
- *6F Guide Kissing*

MediSunik's semi-compliance pre-dilatation catheter which offers numerous features for enhanced crossability, flexibility and extraordinary rewrap to facilitate the physician and procedure. The balloon sizes 1.0 and 1.25mm are designed for small and narrow lesions and for the treatment of CTOs. It is no boundary to cross and no limitation, GO FURTHER AND GO FASTER.

- *Manicured and balanced 0.016" tip*
- *Enhanced balloon flexibility and crossability*
- *5F guiding catheter compatibility and 6F guiding catheter compatibility for kissing balloon technique*
- *Innovative shaft materials with unsurpassed kink resistance*


Balloon Specification

Balloon Diameter(mm)	1.0, 1.25, 1.5, 1.75, 2.0, 2.25, 2.5, 2.75, 3.0, 3.25, 3.5, 3.75, 4.0
Balloon Length(mm)	5, 8, 10, 12, 15, 20, 25, 30
Nominal Pressure	6atm
Rated Burst Pressure	14atm
Average Burst Pressure	20atm
Compliance	Semi-compliance
Balloon Pleats	2/3
Balloon Profile	0.021"($\varnothing$ 1.0mm); 0.030"($\varnothing$ 3.0mm)
Tip Length(mm)	1.5/2.0/2.5
Distal Shaft OD	1.9F
Proximal Shaft OD	2.36F/2.55F/2.7F
Guiding Catheter Compatibility	5F Compatible/6F compatible for kissing
Working Length(cm)	140
Coating	Slider™ Hydro from tip to adjacent guide wire port (external surface); Eel™ coating for guide wire lumen

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

Balloon Dilatation Catheter



- **High Pressure**
- **Non-Compliance**
- **6F Guide Kissing**

MediSunik's truly non-compliance high pressure post-dilatation catheter which offers numerous features for enhanced robust & non-compliance balloon, extraordinary minor balloon axial growth while inflation and 6F guiding catheter compatible. The balloon sizes 2.0-5.0mm and length 8-18mm. The average balloon burst pressure is 30atm. It is satisfied and unlimited, the STRONGER the BETTER and SAFER.

- *Manicured and balanced 0.016" tip*
- *Nominal @ 12atm and Rated Burst Pressure @ 22atm*
- *Minimal Balloon Compliance, 0.55% from nominal to Rated Burst Pressure(RBP)*
- *Axial Balloon Growth, 3% from Nominal to Rated Burst Pressure*
- *Average burst pressure 30atm*
- *5F guiding catheter compatibility and 6F guiding catheter compatibility for kissing balloon technique*
- *Innovative shaft materials with unsurpassed kink resistance*

Balloon Specification

Balloon Diameter(mm)	2.0, 2.25, 2.5, 2.75, 3.0, 3.25, 3.5, 3.75, 4.0, 4.5, 5.0
Balloon Length(mm)	8, 10, 12, 15, 18
Nominal Pressure	12atm
Rated Burst Pressure	22atm(ø2.0-4.0mm) 20atm(ø4.5-5.0mm)
Average Burst Pressure	30atm
Compliance	Non-compliance, 0.55%
Balloon Pleats	3
Tip Length(mm)	2.0/2.5
Distal Shaft OD	2.55/2.6F
Proximal Shaft OD	2.0F
Guiding Catheter Compatibility	5F Compatible/6F compatible for kissing
Working Length(cm)	140
Coating	Slider™ Hydro from tip to adjacent guide wire port (external surface); Eel™ coating for guide wire lumen

Sunny Medical Device (Shenzhen) Co., Ltd.


1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

Disposable Heparin Caps



Ref (Incomplete)	Description
SM-HC-T-MP	Transparent, Medical Paper-plastic Pouch

Disposable towel/sheet Clips



Ref (Incomplete)	Description
SM-TC-W1-MP	White, 1 in 1 Pouch, Medical Paper-plastic Pouch
SM-TC-W2-MP	White, 2 in 1 Pouch, Medical Paper-plastic Pouch
SM-TC-W5-MP	White, 5 in 1 Pouch, Medical Paper-plastic Pouch

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 WhatsApp/Telegram: +852 91216516
Email: sunnymedical_cn@139.com jamesqizhang@gmail.com
Website: www.sunnymedical.com.cn www.sunnymedical.com

CE₂₇₉₇
ISO13485
NMPA



EU Quality Assurance Certificate

Regulation (EU) 2017/745, Annex XI Part A

MDR 761322 R000

Manufacturer: Sunny Medical Device (Shenzhen) Co., Ltd

Address:
1/F and 401 Zhongtianxin Building B
Longgang District
Shenzhen
Guangdong
518172
China

Single Registration Number: CN-MF-000022722

EU Authorised Representative: CMC Medical Devices & Drugs SL

Address:
C/ Morcillo Lengua N18
CP 29006
Málaga
Spain

Scope: See attached Device Schedule

On the basis of our examination of the quality system in accordance with Regulation (EU) 2017/745, Annex XI part A, the quality system meets the requirements of the Regulation. For the placing on the market of Class II and Class IIa devices an additional Annex X certificate is required.

For and on behalf of BSI, a Notified Body for the above Regulation (Notified Body Number 2797):

Graeme Tunbridge
Graeme Tunbridge, Senior Vice President Medical Devices

First Issue Date: 2023-04-27

Current Issue Date: 2023-09-08

Starting Validity Date: 2023-09-08

Expiry Date: 2028-04-26

...making excellence a habit.

Page 1 of 3

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.
This certificate was issued electronically and is bound by the conditions of the contract.

BSI Contact: BSI Group The Netherlands B.V., Sky Building, John H. Kuyperplein 5, 1066 ES Amsterdam, Netherlands. Tel: +31 (0) 20 340 0780
Corporate Contact: BSI Group Assurance Limited, registered in England under number 04325440 at 380 Chiswick High Road, London, W4 4AL, UK.
A Member of the BSI Group of Companies.



By Royal Charter



EU Quality Assurance Certificate

Regulation (EU) 2017/745, Annex XI Part A

MDR 761322 R000

Device Schedule: Class IIa, Custom-made and other devices

Device(s)	Risk Classification
Disposable manifold kits	Class IIa
Manifold and stopcocks	Class IIa
Disposable needles	Class IIa
Pressure lines	Class IIa
Guide Wires	Class IIa
Hemostatic valves	Class IIa
Disposable inflation device kits	Class IIa
Disposable balloon inflation devices	Class IIa
Disposable angioclusion pads	Class IIa
Control syringes	Class IIa
Disposable Thrombore hand/adaptors	Class IIa
Disposable heparin caps	Class IIa
Disposable towel/sheet clips	Class IIa

For Class IIa devices, the Notified Body conformity assessment is limited to the aspects relating to establishing, securing and maintaining sterile conditions.

First Issue Date: 2023-04-27

Current Issue Date: 2023-09-08

Starting Validity Date: 2023-09-08

Expiry Date: 2028-04-26

...making excellence a habit.

Page 2 of 3

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.
This certificate was issued electronically and is bound by the conditions of the contract.

BSI Contact: BSI Group The Netherlands B.V., Sky Building, John H. Kuyperplein 5, 1066 ES Amsterdam, Netherlands. Tel: +31 (0) 20 340 0780
Corporate Contact: BSI Group Assurance Limited, registered in England under number 04325440 at 380 Chiswick High Road, London, W4 4AL, UK.
A Member of the BSI Group of Companies.



EU Quality Assurance Certificate

Regulation (EU) 2017/745, Annex XI Part A

MDR 761322 R000

Certificate History

(References to applicable Common Specifications, Harmonised Standards complied with, and the relevant test and audit reports that support any of the below specifications changes may be requested from Certificates@bsi.com@bsi-group.com)

Date	Reference Number	Action
2023-04-27	3576893	Issued
Current	30007284	Amended = change of EU Authorised Representative to CMC Medical Devices & Drugs SL

First Issue Date: 2023-04-27

Current Issue Date: 2023-09-08

Starting Validity Date: 2023-09-08

Expiry Date: 2028-04-26

...making excellence a habit.

Page 3 of 3

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.
This certificate was issued electronically and is bound by the conditions of the contract.

BSI Contact: BSI Group The Netherlands B.V., Sky Building, John H. Kuyperplein 5, 1066 ES Amsterdam, Netherlands. Tel: +31 (0) 20 340 0780
Corporate Contact: BSI Group Assurance Limited, registered in England under number 04325440 at 380 Chiswick High Road, London, W4 4AL, UK.
A Member of the BSI Group of Companies.



By Royal Charter



Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2016

质量管理体系 ISO 13485:2016

This is to certify that:
兹证明
Sunny Medical Device
(Shenzhen) Co., Ltd.
1/F and 401
Zhongtianxin Building B
Longgang District
Shenzhen
Guangdong
518172
China

深圳市听力医疗设备开发有限公司
914403007488755871
中国
广东省
深圳市
龙岗区
中天信B厂区厂房1楼、401
邮编: 518172

Holds Certificate No: **FM 659381**

持有证书

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 for the following scope:

并运行符合 ISO 13485:2016 要求的质量管理体系, 认证范围如下:

Please see scope page.

For and on behalf of BSI:
BSI 代表:

Graeme Tunbridge
Graeme Tunbridge, Senior Vice President Global Regulatory & Quality

Original Registration Date 首次发证日期: 2016-09-19

Latest Revision Date 最新发证日期: 2024-08-21

Effective Date 生效日期: 2024-10-19

Expiry Date 有效期至: 2027-10-18

Page: 1 of 2



...making excellence a habit.

This certificate was issued electronically and remains the property of BSI and is bound by the conditions of contract. 此证书的所有权属于BSI并受合同条款的限制。

An electronic certificate can be validated at www.bsi.com/bsi-certificates or by telephone +44 (0) 800 730000.
Printed copies can be validated at www.bsi.com/bsi-certificates or by telephone +44 (0) 800 730000.
The information of the certificate can also be found on the official website of the National Certification and Accreditation Administration (www.cnca.gov.cn).

BSI Group The Netherlands B.V. is registered in The Netherlands under number 33343261 | A Member of the BSI Group Holdings B.V.
Information and Contact: BSI Group The Netherlands B.V., John H. Kuyperplein 5, 1066 EP Amsterdam, The Netherlands | Tel: +31 (0) 20 340 780
BSI Group The Netherlands B.V. is registered in The Netherlands under number 33343261 | A Member of the BSI Group Holdings B.V.
(信息在网页及联系方式: BSI Group The Netherlands B.V., John H. Kuyperplein 5, 1066 EP Amsterdam, The Netherlands | 电话: +31 (0) 20 340 780 BSI集团荷兰有限公司, 注册地荷兰, 注册号 33343261 | BSI集团荷兰有限公司在荷兰注册)
China Headquarters (Beijing): 2008 East Ocean Centre, 244 Jiangmenwai, Beijing, 10004, China, Tel: +86 (0) 10 85073000, A Member of the BSI Group of Companies, 中国总部 (北京): 北京中轴线 (中轴线) 东24号东海洋中心2008室 邮编: 100004 电话: +86 (0) 10 85073000 BSI集团中国有限公司。

深圳市听力医疗设备开发有限公司
深圳市龙岗区龙城街道龙平西路4号
中天信B厂区厂房1楼、401
电话： 0755-28957567
邮编： 518172

Sunny Medical Device (Shenzhen) Co., Ltd.
Address: 1/F and 401, Zhongtianxin Building B,
Longgang District, Shenzhen, Guangdong,
518172 China

Tel: +86-755-84532326

WeChat: +86-13823700608

Telegram: +86 138 2370 0608

Website: www.sunnymedical.com.cn

Email: sunnymedical_cn@139.com
jamesqizhang@gmail.com

Sunny Medical Device Limited

Address: FLAT/RM 09A, 05/F, Hong Leong Plaza
33 Lok Yip Road, Fanling, NT Hong Kong

WhatsApp/Telegram: +1 949 365 6488 / +852 91216516

General Manager



Sales Manager



Sunny Catalogue

