


Hawkins™ Blunt Needle Access System

Key Product Attributes

- Provides a lightweight method of targeting, remaining in place while position is confirmed
- The adjustable needle stop localizes the outer cannula to the target, preventing forward migration
- The outer cannula has an echogenic tip, promoting accurate placement under ultrasound guidance
- Numerically ordered centimeter markings on the outer cannula provide reference for depth placement
- Accommodates most guidewires, biopsy needles, and biopsy instruments

Hawkins™ Blunt Needle Access System (Box of 5)

CATALOG NO	DESCRIPTION
216057	16G x 5 cm
216107	16G x 10 cm
216157	16G x 15 cm
218057	18G x 5 cm
218107	18G x 10 cm
218157	18G x 15 cm
218207	18G x 20 cm



Drainage Bags

Drainage Bags are designed for efficient fluid collection during drainage procedures. The 600mL drainage bag has a bottom outlet and an adjustable drain tube length.

Key Product Attributes

- Kink-resistant tube can be cut and fit to the appropriate length
- Flexible tubing allows for patient mobility during drainage
- Twist-type valve provides quick and easy emptying of drainage bag
- Anti-reflux valve prevents retrograde contamination
- Fully adjustable soft belt provides patient convenience

Drainage Bags

CATALOG NO	DESCRIPTION
DBAG600	600mL drainage, with adjustable drain tube length (Box of 20)
DBAG600H	600mL drainage, with adjustable drain tube length (Box of 5)

