

Pro-Mag™ Ultra Automatic Biopsy Instrument

Pro-Mag™ Ultra Automatic Biopsy Instrument is a reusable device that utilizes disposable needles to harvest high quality histological specimens.

Key Product Attributes

- Features a 25mm “throw” advancement
- Trigger is easily actuated (fired) from the front or rear of the instrument
- Built-in safety activates upon cocking to prevent accidental firing
- Easy specimen retrieval without removing the needle from the instrument

Pro-Mag™ Ultra Automatic Biopsy Instrument (Box of 1)

CATALOG NO.	DESCRIPTION
7675	Pro-Mag™ Ultra Automatic Biopsy Instrument

PRO-MAG™ ULTRA



Pro-Mag™ Ultra Biopsy Needles

Specially designed for use with Pro-Mag™ Ultra Automatic Biopsy Instrument. High quality core biopsy needles provide histological specimens for diagnosis of suspect lesions.

Key Product Attributes

- A 19mm sample notch ensures sufficient tissue for clinical diagnosis
- Unique hub design obsoletes needle spacers
- Centimeter markings provide reference for accurate depth placement
- Each needle has an echogenic tip for accurate placement under ultrasound guidance
- Hubs are color-coded for easy gauge size determination
- Optimal co-axial introducer needles available

Pro-Mag™ Ultra Echogenic Biopsy Needles (Box of 10)

CATALOG NO.	DESCRIPTION
765014100	14G x 10cm, optional co-axial needle MCXS1410AX
765014120	14G x 12cm
765014160	14G x 16cm, optional co-axial needle MCXS1416AX
765016100	16G x 10cm, optional co-axial needle MCXS1610AX
765016160	16G x 16cm, optional co-axial needle MCXS1616AX
765016200	16G x 20cm
765018100	18G x 10cm, optional co-axial needle MCXS1810AX
765018160	18G x 16cm, optional co-axial needle MCXS1816AX
765018200	18G x 20cm, optional co-axial needle MCXS1820AX
765018250	18G x 25cm
765018300	18G x 30cm
765020100	20G x 10cm, optional co-axial needle MCXS2010AY
765020160	20G x 16cm, optional co-axial needle MCXS2016AY

Optional Echogenic Co-axial Introducer Needles (Box of 10)

CATALOG NO.	DESCRIPTION
MCXS1410AX	13G x 4.6cm, co-axial to 765014100
MCXS1416AX	13G x 10.6cm, co-axial to 765014160
MCXS1610AX	15G x 4.6cm, co-axial to 765016100
MCXS1616AX	15G x 10.6cm, co-axial to 765016160
MCXS1810AX	17G x 4.6cm, co-axial to 765018100
MCXS1816AX	17G x 10.6cm, co-axial to 765018160
MCXS1820AX	17G x 14.6cm, co-axial to 765018200
MCXS2010AY	19G x 4.9cm, co-axial to 765020100
MCXS2016AY	19G x 10.9cm, co-axial to 765020160

