

Non-Positionable Breast Localization Needles (BLN)

Accura™ BLN

Accura™ Breast Localization Needle (Kopans style) is a traditional, springhook hardwire with stiffener design.

Key Product Attributes

- Skin retention clip prevents forward migration of wire during patient transport
- Wire is marked in two places: I indicates wire at distal tip of needle; II indicates barb deployment
- End hole design allows for dye injection or aspiration of fluid

Accura™ Breast Localization Needles (Kopans style) (Box of 10)

CATALOG NO.	DESCRIPTION
BLN2003	20G x 3 cm
BLN2005	20G x 5 cm
BLN20075	20G x 7.5 cm
BLN2010	20G x 10 cm
BLN2015	20G x 15 cm
BLN2103	21G x 3 cm
BLN2105	20G x 5 cm
BLN21075	20G x 7.5 cm
BLN2110	20G x 10 cm



Accura™ II BLN

Accura™ II Breast Localization Needle (Sadowsky Style) comes with a sterile separately packaged 23G stiffening cannula. Easily placed over the wire by the surgeon, the cannula indicates the course of the wire, provides a palpable landmark, and prevents accidental wire transection.

Key Product Attributes

- Centimeter markings provide reference for depth placement
- Throated needle hub for easy wire loading
- Wire is marked in two places: I indicates wire at distal tip of needle; II indicates barb deployment

Accura™ II Breast Localization Needles (Sadowsky style) (Box of 10)

CATALOG NO.	DESCRIPTION
266050	20G x 5 cm with 23G stiffening cannula
266075	20G x 7.5 cm with 23G stiffening cannula
266100	20G x 10 cm with 23G stiffening cannula
266125	20G x 12.5 cm with 23G stiffening cannula

