

# Bone & Bone Marrow Biopsy

## T-Lok™ Bone Marrow Biopsy Needles

T-Lok™ Bone Marrow Biopsy Needles are designed to provide quality biopsies by entrapping bone marrow in the needle without altering the sample.

### Key Product Attributes

- Unique cannula design provides outstanding quality specimens
- Ergonomically designed twist-lock handle provides greater clinician comfort
- Double-diamond stylet point easily penetrates hard bone
- Innovative design provides an exceptional ability to capture and retain bone marrow specimens
- Eliminates painful deflection technique to harvest a diagnostic quality specimen

T-Lok™ Bone Marrow Biopsy Needles (Box of 10)	
CATALOG NO.	DESCRIPTION
DBMNJ0804TL	8G x 4" with extraction cannula
DBMNJ0806TL	8G x 6" with extraction cannula
DBMNJ1104TL	11G x 4" with extraction cannula
DBMNJ1106TL	11G x 6" with extraction cannula
DBMNJ13025	13G x 2.5" - no extraction cannula
DBMNJ1304	13G x 4" - no extraction cannula



## Bone Marrow Aspiration Needles

Safe aspirations at the sternum can be performed by adjusting the needle stop to control the depth of penetration. The needle stop can be removed if aspirating at the iliac crest.

### Key Product Attributes

- Luer lock connector on the handle for secure syringe attachment
- An adjustable, removable needle stop safely controls depth of penetration

Bone Marrow Aspiration Needles (Box of 10)	
CATALOG NO.	DESCRIPTION
DBMNI1501	15G x 2.688" max.
DBMNI1504	15G x 4" max.
DBMNI1601	16G x 2.688" max.

