

## Repositionable Breast Localization Needles (BLN) (cont.)

### Hawkins™ II BLN

Hawkins™ II BLN is repositionable to achieve accurate placement.

#### Key Product Attributes

- Available in two wire choices: FlexStrand™, a 19-strand cable for maximum flexibility with minimal risk of transection, or a traditional hardwire, which minimizes flex during placement
- Both wires offer a side barb for temporary fixation during a confirming mammogram. The needle is easily repositioned, if necessary
- Easy-to-read centimeter markings facilitate depth placement
- Throated needle hub for easy wire loading
- Skin retention clip prevents forward migration of wire during patient transport
- Wire is marked in three places: I indicates wire at distal tip of needle; II indicates side barb deployment; and III indicates distal barb deployment
- Standard marking cannula on wire, beginning 2.5cm from distal tip, serves as a landmark for the surgeon

#### Hawkins™ II Breast Localization Needles w/FlexStrand™ (Box of 10)

CATALOG NO.	DESCRIPTION
252050	20G x 5cm
252075	20G x 7.5cm
252100	20G x 10cm
252125	20G x 12.5cm

#### Hawkins™ II Breast Localization Needles w/Hardwire (Box of 10)

CATALOG NO.	DESCRIPTION
242050	20G x 5cm
242075	20G x 7.5cm
242100	20G x 10cm
242125	20G x 12.5cm

