

Guidewires

Guidewires - Fixed Core Straight

Stainless Steel			
CATALOG NO.	DIAMETER	LENGTH	UNITS
395451	.018"	40cm	10
395200	.018"	80cm	10
395201	.018"	100cm	10
395536	.021"	25cm	10
395513	.021"	40cm	10
395203	.021"	80cm	10
395205	.021"	125cm	10
395509	.021"	150cm	10
395206	.025"	40cm	10
395711	.035"	75cm	10

PTFE Coated Stainless Steel			
CATALOG NO.	DIAMETER	LENGTH	UNITS
395225	.025"	150cm	10
395687	.035"	150cm	10
395234	.038"	150cm	10

Exchange Wires - PTFE Coated			
CATALOG NO.	DIAMETER	LENGTH	UNITS
395367	.025"	260cm	10
395369	.035"	260cm	10
395370	.038"	260cm	10

Mullins Tip Deflector Wires - Stainless Steel			
CATALOG NO.	DIAMETER	LENGTH	UNITS
395945	.018"	100cm	10
395948	.021"	100cm	10

Guidewires - Fixed Core J Tip

PTFE Coated Stainless Steel					
CATALOG NO.	DIAMETER	LENGTH	TIP CURVE RADIUS	SPECIAL FEATURE	UNITS
395420	.025"	150cm	1.5mm		10
395787	.035"	150cm	1.5mm		10
395523	.021"	150cm	3mm		10
395559	.025"	150cm	3mm		10
395600	.032"	150cm	3mm		10
395282	.035"	150cm	3mm		10
388282	.035"	150cm	3mm	Classic	10
393282	.035"	150cm	3mm	Flexible	10
393182	.035"	150cm	3mm	Z-Tip	10
390282	.035"	150cm	3mm	XTB	10
393186	.035"	175cm	3mm	Z-Tip	10
390186	.035"	175cm	3mm	XTB	10
395284	.038"	150cm	3mm		10
388284	.038"	150cm	3mm	Classic	10
393184	.038"	150cm	3mm	Z-Tip	10
390284	.038"	150cm	3mm	XTB	10
393187	.038"	175cm	3mm	Z-Tip	10

PTFE Coated Stainless Steel Double Ended				
CATALOG NO.	DIAMETER	LENGTH	TIP CURVE RADIUS	UNITS
395998	.035"	50cm	Str & 3mm	10
395919	.035"	150cm	Str & 1.5mm	10

Exchange Wires - PTFE Coated				
CATALOG NO.	DIAMETER	LENGTH	TIP CURVE RADIUS	UNITS
395585	.025"	260cm	3mm	10
395722	.035"	200cm	3mm	10
395794	.035"	260cm	3mm	10

Guidewires - Specialty

Stainless Steel				
CATALOG NO.	DIAMETER	LENGTH	TIP CURVE RADIUS	UNITS
395510	.021"	40cm	3mm	10
395573	.025"	40cm	3mm	10
395271	.035"	40cm	3mm	10
395654	.035"	50cm	3mm	10
395752	.035"	75cm	3mm	10

Bentson Wires - PTFE Coated				
CATALOG NO.	DIAMETER	LENGTH	SPECIAL FEATURE	UNITS
395930E	.035"	150cm		10
395170	.035"	180cm		10
390033	.035"	260cm		10
395933C	.038"	150cm		10

Stainless Steel Double Ended				
CATALOG NO.	DIAMETER	LENGTH	TIP CURVE RADIUS	UNITS
395415	.025"	65cm	Str & 3mm	10
395539	.032"	65cm	Str & 3mm	10
395680	.035"	45cm	Str & 3mm	10
395899	.035"	65cm	Str & 3mm	10
395868	.035"	80cm	Str & 3mm	10

Movable Core PTFE Coated					
CATALOG NO.	DIAMETER	LENGTH	TIP CURVE RADIUS	SPECIAL FEATURE	UNITS
395258	.035"	150cm	Straight		10
395309	.032"	150cm	3mm		10
395721	.035"	145cm	3mm	Tapered	10
395320	.035"	150cm	3mm		10
388300	.035"	150cm	3mm	Tapered	10
395312	.035"	150cm	3mm	Classic	10
395315	.038"	150cm	3mm		10

.018" Access Wires

Guidewires come in two lengths and are available with stainless steel or nitinol cores. All guidewires have a flexible radiopaque tip.

Key Product Attributes

- Smooth transition from rigid core to radiopaque coil tip provides excellent maneuverability
- Nitinol offers kink-resistance and torque control
- Stainless steel for maximum back-up

Guidewire (Box of 5)	
CATALOG NO.	DESCRIPTION
107018080	.018" X 80cm stainless steel guidewire
190118040	.018" X 40cm nitinol guidewire
190118080	.018" X 80cm nitinol guidewire