



# WORKER<sup>®</sup> AMPLATZ

GUIDEWIRE




# WORKER AMPLATZ GUIDEWIRE

The Right Combination of Stiffness, Durability, and Rail Support



- **Large Stainless Steel (SS) Core with SS Flat-Wire Coils** increases rigidity, tactile feedback, and support while advancing catheters through the vasculature
- **PTFE Coating** reduces friction and increases lubricity along the entire surface
- **Gradual Taper with Core-to-Tip Support and Flexibility** for atraumatic navigation through tortuous anatomy

## Superior Stiffness Versus the Competition\*



\*Data on File

## Ordering Information

Worker Amplatz Guidewire - Straight			
Item #	Diameter	Wire Length	Tip Length
124135080	0.035"	80cm	3.5cm
124135150	0.035"	150cm	3.5cm
124135180	0.035"	180cm	3.5cm
124135260	0.035"	260cm	3.5cm
124035090	0.035"	90cm	7.5cm
124035150	0.035"	150cm	7.5cm
124035180	0.035"	180cm	7.5cm
124035260	0.035"	260cm	7.5cm

Worker Amplatz Guidewire - J-Curve			
Item #	Diameter	Wire Length	Tip Length
124235090	0.035"	90cm	7.5cm
124235150	0.035"	150cm	7.5cm

Worker is a registered trademark of Argon Medical Devices, Inc.

2600 Dallas Parkway Ste. 440, Frisco, TX 75034  
 P: 800.927.4669 • F: 903.677.9396 • argonmedical.com



96-6400-49/0419A  
 Copyright 2019 Argon Medical Devices, Inc.